

# Abdominal Binder, 10"

With Tapered Design

Model 8619

# Grafco®

**Manufactured for:**  
Graham-Field Health Products  
Atlanta, Georgia 30360  
Made In U.S.A.  
[www.grahamfield.com](http://www.grahamfield.com)

## Indications

Ideal for providing compression and support for strains and weakness of the abdominal area. Provides abdominal and lumbar support. Abdominal, post natal or post surgical application. Easy adjustment with Velcro closure.

## Features and Benefits

- 9" one-panel elastic with plush inner body
- Tapered design for better fit
- Effortless Velcro hook and loop
- Options: Available in 6", 9", 10", or 12" height
- Latex Free



## Instructions:



First, open the binder to its full length. With the Velcro fastener on the patient's right and the two seams centered at the patient's back, position the binder under the patient with the soft lining against the skin. Position the lower edge crossing halfway down the buttocks and is low on the patient's pelvic area.



From the patient's left side, stretch out the end of the binder with one hand and position it across the stomach, holding it firmly in position.



With the other hand, grasp the other end of the binder in the middle, stretching it over the patient's body and firmly engage the Velcro surfaces.



To adjust the fit, simply peel the surfaces apart and reposition them. The closure can be positioned to conform to any patient's body shape.

## Please Note:

An abdominal binder should be used in conjunction with a medical professional's direction. We cannot guarantee wearing an abdominal binder will promote healing.

## Warning:

If you experience any pain, swelling, sensation changes, or any unusual reactions while using this product, consult your medical professional immediately.

*The abdominal binder may be machine or hand washed in a mild detergent. Warm or cold water temperature is recommended. Before washing make sure that the Velcro surfaces are securely fastened to prevent lint accumulation which can interfere with the Velcro fibers and weaken their locking properties. The binder may be hung dry.*

## Sizing:

To insure you select the proper size that will fit you best, select the size that represents the measurement of the fullest part of the abdomen. In most cases, this will insure the most proper fit. However, when the waist size and hip size differ 8" or more, select a size that fits in the middle of the range between the waist size and hip size.

### Size Chart

Large  
XL  
2XL

### Fits

36" - 42"  
42" - 48"  
48" - 56"

**Proudly Made  
in the USA!**



LATEX FREE PRODUCT