





Note: The most current version of these instructions can be found online at www.grahamfield.com

- ⚠ **WARNING: Important! Read and understand these instructions before assembling or using the Padded Hammock Sling. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Padded Hammock Sling is not properly applied, personal injury and damage to the Padded Hammock Sling could result.**
- ⚠ **WARNING: If components are damaged or missing, contact your dealer immediately.**
- ⚠ **WARNING: Carefully inspect the sling weekly for wear and damage to seams, fabric, straps, and strap loops. Torn, cut, frayed or broken slings can fail, resulting in serious personal injury to the user. Use only slings that are in good condition. Discard and destroy old, unusable slings.**
- ⚠ **WARNING: The Padded Hammock Sling has a maximum weight capacity of 400 lb. DO NOT exceed the maximum weight capacity.**
- ⚠ **WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper application or use of this product.**

APPLICATION

1. Check the patient's weight. Ensure that the patient's weight does not exceed the Hammock Sling's maximum weight capacity.
2. Review the following procedures carefully before attempting lift; call Graham-Field Technical Support at 800-347-5678 with any questions.
3. See the Sling Loop Options table at right to determine the appropriate position.
4. Position the sling under and behind the patient's back, with color adjustment strap on outside, and straps **A** and **B** at shoulder height.
5. Pull leg straps **E** and **F** down behind patient and between patient's legs, as shown in Figure 2.
6. Hook leg straps **E** and **F** through small loop **G**, as shown in Figure 3.
7. Slip strap ends **C** and **D** through strap loop ends **E** and **F** and hook to spreader bar at **e** and **f** as shown in Figures 2 and 4.
8. Do not lock caster brakes, in order to allow the lift to properly balance weight.
9. Proceed with lift and transfer as shown in Figure 5.

Sling Loop Options	
Long loop	reclined position
Center loop	semi-reclined position
Short loop	most upright position



MAINTENANCE

To clean Hammock Sling, use a mild soap and water solution and a clean cloth. Rinse with clear water, then allow to air dry.

▲ **Caution: Never use solvents, abrasive cleaners, or scouring pads on any part of the Hammock Sling.**

WARRANTY

GF Health Products, Inc. offers a one-year limited warranty against manufacturer's defects on the Padded Hammock Sling. If a product is deemed to be under warranty, GF Health Products, Inc. shall provide, at its option, (1) replacement of any defective part or product or (2) a credit of the original selling price made to GF Health Products, Inc.'s initial customer. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.


FEATURES

Padded Hammock Sling

- Ideal for single or double amputees and patients with limited muscle control
- Provides head, neck, and full body support
- Designed to eliminate pressure points found in other sling styles
- Deluxe padding for maximum comfort
- 400 lb maximum weight capacity

HP222 Medium 1 ea

HP223 Large 1 ea



ITEM	Sling Size	Body Width	Body Length	Total Width	Total Length
HP222	Medium	30 1/2"	21"	74 3/4"	39 1/2"
HP223	Large	33 1/4"	24 1/4"	76"	47 1/4"



GRAHAM-FIELD

GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice. Graham-Field and Lumex are registered trademarks of GF Health Products, Inc.

© October 2009 GF Health Products, Inc. • 2935 Northeast Parkway • Atlanta, GA 30360 • tel 800-347-5678 • fax 800-721-0601 • www.grahamfield.com