

*Note: The most current version of these instructions can be found online at [www.grahamfield.com](http://www.grahamfield.com)*


- ⚠ **WARNING: Important! Read and understand these instructions before assembling or using the Full-Body Sling. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Full-Body Sling is not properly applied, personal injury and damage to the Full-Body Sling could result.**
- ⚠ **WARNING: If components are damaged or missing, contact your dealer immediately.**
- ⚠ **WARNING: Carefully inspect the sling weekly for wear and damage to seams, fabric, straps, and strap loops. Torn, cut, frayed or broken slings can fail, resulting in serious personal injury to the user. Use only slings that are in good condition. Discard and destroy old, unusable slings.**
- ⚠ **WARNING: DO NOT exceed the Full-Body Sling's maximum weight capacity.**
- ⚠ **WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper application or use of this product.**

### APPLICATION

1. Check the patient's weight and the sling's maximum weight capacity. Ensure that the patient's weight does not exceed the sling's maximum weight capacity.
2. Review the following procedures carefully before attempting lift; call Graham-Field Technical Support at 800-347-5678 with any questions.
3. See the Sling Loop Options table at right to determine the appropriate position.
4. See Figure 1 for correct positioning of sling on spreader bar.
5. Position the sling under and behind the patient's back with the color adjustment strap on outside. Drape **A** and **B** over the shoulders as shown in Figure 2.
6. Carefully, with help from another caregiver, if needed, pull the sling portion **C** and **D** forward to the bend of the knee as shown in Figure 2.
7. Hook sling loops **A** and **B** to spreader bar **a** and **b** as shown in Figure 3.
8. Hook sling loops **C** and **D** to spreader bar **c** and **d** as shown in Figure 3.
9. Do not lock caster brakes, in order to allow the lift to properly balance weight.



Sling Loop Options	
Long loop	reclined position
Center loop	semi-reclined position
Short loop	most upright position



10. Proceed with lift and transfer as shown in Figure 4.

**MAINTENANCE**

To clean sling, use a mild soap and water solution and a clean cloth. Rinse with clear water, then allow to air dry.

▲ **Caution: Never use solvents, abrasive cleaners, or scouring pads on any part of the sling.**

**WARRANTY**

GF Health Products, Inc. offers a one-year limited warranty against manufacturer's defects on the Full-Body Sling. If a product is deemed to be under warranty, GF Health Products, Inc. shall provide, at its option, (1) replacement of any defective part or product or (2) a credit of the original selling price made to GF Health Products, Inc.'s initial customer. The warranty does not include any labor charges incurred in replacement part(s) installation or any associated freight or shipping charges to GF Health Products, Inc.



**FEATURES**

**Mesh Full-Body Commode Sling**

Quick-draining/drying mesh fabric offers full head and neck support. For use with 4-point spreader option. Designed to meet the requirements of HCPCS Code E0621.

- DSLRL114 Medium 1 ea
- DSLRL115 Large 1 ea
- DSLRL116 X Large 1 ea
- DSLRL141 Heavy Duty, XX Large 1 ea

**Mesh Full-Body Sling**




Quick-draining/drying mesh fabric offers full head and neck support. For use with 4-point spreader option. Designed to meet the requirements of HCPCS Code E0621.

- DSLRL110 Medium 1 ea
- DSLRL111 Large 1 ea
- DSLRL140 Heavy Duty, XX Large 1 ea

**Solid Full-Body Sling**

Fully-padded fabric offers a soft, yet support surface. Offers full head and neck support. For use with 4-point spreader option. Designed to meet the requirements of HCPCS Code E0621.

- DSLRL112 Medium 1 ea
- DSLRL113 Large 1 ea
- DSLRL117 X Large 1 ea



ITEM	Size	Body Width	Body Length	Commode Opening to Top of Back (DSLRL114, DSLRL115, DSLRL116, DSLRL141 only)	Maximum Weight Capacity
DSLRL110, DSLRL112, DSLRL114	Medium	42 3/4"	50 3/4"	34 3/4"	450 lb
DSLRL111, DSLRL113, DSLRL115	Large	44 3/4"	54 1/4"	38"	450 lb
DSLRL116, DSLRL117	X Large	44 3/4"	60 1/4"	43"	450 lb
DSLRL140, DSLRL141	Heavy Duty, XX Large	42"	51"	41 1/2"	600 lb



**GRAHAM-FIELD**

GF Health Products, Inc. is not responsible for typographical errors. Packaging, warranties, products and specifications are subject to change without notice. Graham-Field and Lumex are registered trademarks of GF Health Products, Inc.

© October 2009 GF Health Products, Inc. • 2935 Northeast Parkway • Atlanta, GA 30360 • tel 800-347-5678 • fax 800-721-0601 • www.grahamfield.com