



Figure 1



Figure 2

### Tools Required:

Phillips Screwdriver

### Table Hanger Installation:

1. Lock rear casters.
2. Place table hanger assembly over back cross brace. (See Figure 1)
3. Insert carriage bolts through square hole in rear of hanger bracket. Add washers and wing nuts. Tighten wing nuts.

### Support Strap Installation:

1. Remove upholstery from back of chair.
2. Insert plastic rivets through grommets on ends of strap into small holes in chair back bend.
3. Push or tap lightly with small hammer on plunger of rivet until flush, this will lock rivet in place. (See Figure 2)





Figure 3



Figure 4



Figure 5



Figure 6

### Receiving Tubes Installation:

1. Standing behind the chair, tip it over on its left side.
2. Remove two front screws and washers from underneath the right arm. With the tube facing front, reinsert washer and screws through the plate, back into the armrests (See Figure 3).
3. Tip chair onto other side. Repeat previous step with long tube. Pull pull pin and slide lockout forward (See figure 3).
4. Put chair upright.
5. Align receiving tubes to table by holding table at handles and rotate round tubes until they can be inserted into square tubes. Holding handles, slide table in and out a few times to check alignment (See Figure 4).

Push on table center while tightening screws under armrests.

6. If table still feels stiff, loosen screws by two turns and readjust plates.

### Table Lock:

To make the table lock, pull the pull button down and slide the lockout back (See Figure 5).

### To Store Table:

1. Remove tray from receiving tubes and push handles underneath tray table.
2. Place edge of the table on hanger.
3. Slide strap over top of table (See Figure 6).



**Graham-Field Health Products**

2935 Northeast Parkway, Atlanta, Georgia 30360

telephone: 800-347-5678, 770-447-1609

fax: 800-726-0601, 678-291-3232

<http://www.grahamfield.com>