



PLEASE SAVE THESE INSTRUCTIONS FOR FUTURE USE

Tools required: #2 Phillips Screwdriver

Info: The most current version of these instructions can be found online at www.grahamfield.com.

- ⚠ WARNING: Important! Read and understand these instructions before assembling or using the Lumex Everyday Hip Chair with Adjustable Footrest. If you do not understand any part of these warnings or instructions, contact a healthcare professional for direction in the use of this product. If the Hip Chair is not properly assembled, personal injury and damage to the Hip Chair could result.**
- ⚠ WARNING: If components are damaged or missing, contact your GF distributor immediately. DO NOT use substitute parts. Use only Lumex replacement parts. Non-Lumex replacement parts could cause personal injury and damage to the Hip Chair.**
- ⚠ WARNING: This product has a maximum weight capacity of 350 lb, EVENLY DISTRIBUTED.**
- ⚠ WARNING: Notice for California Customers- California Proposition 65
WARNING: This product contains a chemical known to the State of California to cause cancer and reproductive or developmental harm.**
- ⚠ WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper assembly, installation or use of this product.**

INTENDED USE

The Everyday Hip Chair with Adjustable Footrest (shown on next page) is intended for use by those who have difficulty sitting down or standing up, including those with arthritis or who have undergone hip or knee surgery.

CE CERTIFICATION

 **GF Health Products, Inc.**
2935 Northeast Parkway Atlanta, Georgia 30360
Tel: 678-291-3207 Web: www.grahamfield.com

 EC Authorized Representative
CEpartner4U BV, ESDOORNLAAN 13
3951 DB MAARN The Netherlands
+31 (0)6-516.536.26
www.cepartner4u.com



ASSEMBLY (REFER TO PICTURE AT RIGHT)

▲ NOTICE: When installing screws during assembly, always start the screws by hand, then tighten them with the screwdriver. Failure to start the screws by hand could cause damage to the Hip Chair.


1. Remove the frame and seat-back assembly from the carton. Stand the Hip Chair upright on all four legs.
2. Taking care to support the seat-back assembly so it does not fall when the screws are removed, remove the four screws that secure the seat to the seat support underside.
3. Position the seat-back assembly on top of the seat support with the seat facing forward toward the footrest as shown at right.
4. **Using the long screws in the seat front and the short screws in the seat rear**, reinstall the four screws through the frame underside into the seat underside. Tighten the screws.
5. Screw the footrest height adjustment knobs into the threaded bushings in the front legs until firmly seated — permit the thread locking compound on the knobs to cure at least 24 hours for maximum retention before using the chair.
6. Position the footrest as shown at right and insert its posts into the front leg sockets while operating the height adjustment knobs as described in the FOOTREST HEIGHT ADJUSTMENT steps that follow.



Hip Chair with Adjustable Footrest, Shown Assembled

FOOTREST HEIGHT ADJUSTMENT

The Lumex Everyday Hip Chair has an adjustable-height footrest to make the chair more comfortable depending on the patient's height requirements. Refer to pictures that follow to adjust the footrest.



1. Pull either footrest height adjustment knob outward and rotate it one quarter-turn to lock the knob in the outward position.
2. Pull the opposite knob outward and rotate it one quarter-turn to lock the knob in the outward position.
3. Adjust the footrest up or down to the desired height.
4. Once the desired footrest height is reached, turn either knob back to unlock the knob and lock the footrest at the desired height.
5. Turn the opposite knob to unlock it and lock the footrest at the desired height.

6. Move the footrest up or down to ensure the footrest is locked in desired position; if not, repeat steps.
- ⚠ **WARNING: DO NOT stand on the footrest as the chair may tip. The footrest is designed to be used as a rest for the legs, but not designed to support the total body weight of an individual.**
- ⚠ **WARNING: DO NOT use the chair until all assembly steps have been completed as described and footrest is locked in place at an appropriate height for the user.**

MAINTENANCE AND CLEANING

To clean Hip Chair, wash with a clean, soft cloth dampened with soap and water or a mild household detergent, wipe with a clean, soft cloth dampened with clear water, and dry with a clean, soft cloth.

⚠ **WARNING: Inspect Hip Chair weekly; tighten fasteners as needed to ensure secure installation.**

▲ **NOTICE: Never use abrasive cleaners or scouring pads on any part of the Hip Chair.**

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. (“GF”) warrants to the original purchaser only that it will replace or repair components, at GF’s sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported to GF’s customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable.

The warranted components and time period are set forth below:

Frame: two years
Upholstery: one year

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

A GF Customer Service Representative must authorize warranty service. Please contact the GF Customer Service department by calling 678-291-3207, sending a fax request to 770-368-2386 or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;

- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: casters, filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 6) Any labor or shipping charges incurred in the replacement part installation or repair;
- 7) Costs and expenses of regular maintenance and cleaning; and
- 8) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS.

THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply.
- 2) Freight claims must be notated on the Bill of Lading and must be made with immediacy. The ICC regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within thirty (30) days of the invoice date.



GRAHAM-FIELD
www.grahamfield.com

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this literature are based on the latest product information available at the time of printing. The most current product information, including the most current version of these instructions, can be found online at www.grahamfield.com. **Graham-Field** and **Lumex** are registered trademarks of GF Health Products, Inc.

© 2012 GF Health Products, Inc. • 2935 Northeast Parkway, Atlanta, GA 30360 • tel 770-368-4700 • Made in China