Manufactured for GF Health Products, Inc.
 Atlanta GA 30360
 770-368-4700

GF Health Products, Inc. is not responsible for typographical errors. All illustrations, specifications, packaging and warranties contained in this literature are based on the latest product information available at the time of printing. The most current product information, including the most current version of this manual, can be found online at www.grahamfield.com.


Graham-Field, Lumiscope, and Lumiscope For The Quality Of Life are registered trademarks of GF Health Products, Inc. © 2014 GF Health Products, Inc.

Model 1146
Digital Blood Pressure Monitor
Fully Automatic Wrist Style

Instruction Manual *

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1146-INS-LAB-RevA14

Thank you for choosing the Lumiscope® Digital Blood Pressure Monitor!
 Please read this instruction manual before using your Blood Pressure Monitor.
 Please keep this instruction manual safe for future use.

Table of Contents

Before Using the Monitor

- Introduction 1
- Symbols 2
- Important Safety Precautions 3-4
- Know Your Blood Pressure Monitor 5-6
- IHB Introduction 7
- Before Taking a Measurement 8

Operating Instructions

- Battery Installation / Replacement 9
- Adjusting Settings 10
 - 1. To Set Year 10
 - 2. To Set Date and Time 11-12
 - 3. To Set positioning sensor 12
 - 4. To Set Display of the Date and Time 13-14
 - 5. To Set Volume 15
 - 6. To Set Language (English / Spanish) 17
- Proper Measurement Method 18
- Take a Measurement 19-21

Table of Contents

- Use the Memory Function 22
 - To read the average value 22
 - To read the measurement value 23
 - To delete the memory value 24
- About Blood Pressure 25
 - Blood Circulation 25
 - Health and Blood Pressure 26
 - Classification of Blood Pressure 27-28

Care and Maintenance

- Troubleshooting 29
 - Error Indicators 29
 - Eliminate Fault 30
- Care and Maintenance 31
 - Care 31
 - Maintenance 32
- Specifications 33-36

Warranty

- Limited Warranty 37-38
- Appendix A: EMC (Electro-Magnetic-Compatibility)
 - EMC Information 39-44

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Introduction

Thank you for purchasing the Lumiscope Model 1146 Fully Automatic Wrist Style Digital Blood Pressure Monitor.

Intended Use: This device is intended to measure blood pressure and pulse rate of adults only, at home or medical center. The device is NOT recommended for neonatal use.

The unit uses the oscillometric method of blood pressure measurement, which means the unit detects the movement of your blood through your brachial artery, and converts your blood pressure into a digital reading. The unit is simple to use because a stethoscope is not needed while using an oscillometric monitor.

The unit automatically stores 180 (2 * 90) sets of measurement values with the measuring date and time. You can read the stored data conveniently by pressing the memory button.

The unit comes with the following components:

- Main Unit
- Storage Case
- Instruction Manual printed in English and Spanish

The unit has the following advanced functions:

- Backlight
- Voice Function
- Positioning Sensor
- Display of the date and time

1

Symbols

SYMBOL KEY	
 WARNING	Indicates a potentially hazardous situation which, if not avoided, could result in death or serious personal injury.
 CAUTION	Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate personal injury to the user or patient.
 NOTICE	Indicates a potentially hazardous situation which, if not avoided, could result in damage to the equipment or other property.
	Type B applied part
	Refer to instruction manual/booklet
	SERIAL NUMBER
	The marking of electrical and electronics devices according to Directive 2002/96/EC. The device, accessories and packaging must be disposed of correctly at the end of their usage. Please follow Local Ordinances and Regulations for disposal.

2

Important Safety Precautions

SAFETY GUIDELINES - PLEASE READ BEFORE USE

-  **Important!** Read and understand these instructions before using the Blood Pressure Monitor. If you do not understand any part of these instructions, contact your medical professional or Graham-Field dealer for direction in the use of this product.
-  **WARNING:** If components are damaged or missing, contact your Graham-Field dealer immediately. DO NOT use substitute parts. Use only Lumiscope replacement parts. Non-Lumiscope replacement parts could cause personal injury and damage to the Blood Pressure Monitor.
-  **WARNING:** Contact your physician for specific information about your blood pressure. Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your healthcare provider.
-  **WARNING:** Operate the device only as intended. DO NOT use it for any other purpose.
-  **WARNING:** Do not use the device where flammable gas (such as anesthetic gas, oxygen or hydrogen) or flammable liquid (such as alcohol) are present.
-  **WARNING:** Notice for California Customers- California Proposition 65 **WARNING:** This product contains a chemical known to the State of California to cause cancer and reproductive or developmental harm.
-  **WARNING:** GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation or use of this product.


3

Important Safety Precautions

-  **NOTICE:** Do not use a cellular phone near the device. It may result in operational failure.
-  **NOTICE:** Avoid use of the device in areas of high heat; this may affect the accuracy of the measurement.
-  **NOTICE:** Do not modify, disassemble, or attempt to repair the device or components.


4

Know Your Blood Pressure Monitor

5

Know Your Blood Pressure Monitor


6

IHB Introduction

Irregular Heartbeat Symbol (IHB)

When the unit detects an irregular rhythm two or more times during the measurement, the Irregular Heartbeat Symbol will appear on the display with the measurement values.

An irregular heartbeat rhythm is defined as a rhythm that varies by less than 25% from the average rhythm or 25% more than the average rhythm detected while the monitor is measuring the systolic and diastolic blood pressure.

If the Irregular Heartbeat Symbol (IHB) displays with your measurement results, we recommend you consult your physician. Follow the directions of your doctor.

7

Before Taking a Measurement

1. Avoid eating, drinking alcohol, smoking, exercising and bathing for 30 minutes before taking a measurement.
2. Stress raises blood pressure. Avoid taking measurements during stressful times.
3. The cuff can be applied to your left or right wrist. Always measure your blood pressure on the same wrist.
4. Measurements should be taken in a quiet place.
5. Position the unit at heart level throughout the measurement.
6. Remain still and do not talk during the measurement.
7. Keep a record of your blood pressure and pulse readings for your physician. A single measurement does not provide an accurate indication of your true blood pressure. You need to take and record several readings over a period of time. Try to measure your blood pressure at the same time each day for consistency.
8. Wait 30 ~ 60 seconds between measurements. The wait time allows the arteries to return to their condition prior to taking the blood pressure measurement. You may need to increase the wait time depending on your individual physiological characteristics.

8

Battery Installation / Replacement

- Remove the battery cover from the battery compartment.
Use two identical 1.5V alkaline batteries, type AAA / LR03.
- Install the batteries so the + (positive) and - (negative) battery polarities match the polarities of the battery compartment as indicated.

- Replace the battery cover.


- ▲ NOTICE: Use two identical 1.5V alkaline batteries type AAA / LR03.
- ▲ NOTICE: DO NOT use rechargeable batteries.
- ▲ NOTICE: Replace BOTH batteries if the low battery symbol displays.
- ▲ NOTICE: DO NOT leave low batteries in the battery compartment! They may leak and cause damage to the unit.

9

Adjusting Settings

▶ 1 To set year

- All information will appear on the display when the batteries are installed for the first time. After 2 seconds the year flashes on the display. Press the "SET" button and continue with step 2.
- The year appears and flashes after pressing the "SET" button for 3~4 seconds while the unit is off. Press the "SET" button again and continue with step 2.
Press the "MEM" button one time to increase one year, twice for two years, etc. Keep the "MEM" button depressed and the value quickly increases.

The year can be set between 2012 and 2030.
After displaying 2030, the year will return to 2009.


10

Adjusting Settings


▶ 2 To set date and time

Press the "MEM" button to change the month, then press the "SET" button and continue with next step. Then set the day, the time system conversion, the hour and the minute by using the same method. Continue with step 3.


Month Setting


Day Setting

time system conversion


Hour Setting

11


Adjusting Settings

Minute Setting

▶ 3 To set positioning sensor

The Positioning Sensor symbol "⌂" will appear. Press the "MEM" button to change the current setting. If the positioning system is turned on, the LCD will display "on". If the positioning system is turned off, the LCD will display "off". Press the "SET" button and continue with step 4.


12

Adjusting Settings

► 4 To set Display of date and time

Press the "MEM" button and confirm whether the unit displays the date and time after turning off.

Press the "SET" button and continue with step 5.


When the display function of the date and time is ON, the LCD will display "on".
When the display function of the date and time is OFF, the LCD will display "off".

13

Adjusting Settings

If selected to display the date and time, the current date and time will display when unit has been turned off.


14

Adjusting Settings

► 5 To set Volume

Press the "MEM" button and the volume will cycle from 0 to 7.
Select desired volume level, press the "SET" button, and continue to step 6.
Press the "START/STOP" button to save the current setting and turn the unit off.

When the volume is level 0, the LCD will display "off" and the voice function is OFF.

When the volume is set between level 1 and level 7, the LCD will display "on" and the voice function is ON.

15

Adjusting Settings

► 6 To set Language (English / Spanish)

After Volume has been set, "0 mmHg" will display. Press the "SET" button again to continue to the Set Language screen.

Press the "MEM" button to cycle between languages (English and Spanish). Choose "ESP" to select SPANISH: the talking language will change to Spanish.

Note that the language setting can only be changed during the initial settings adjustment. The default language is English. If you later wish to change the language setting, the batteries must be removed and steps 1 to 6 repeated.

After Language has been set, "0 mmHg" will display.
Press the "SET" button again to continue to the set user screen.


16

Adjusting Settings

7 To set user

Press the "SET" button and the user symbol will appear when turned off or the display function of the date and time is ON. Then press the "MEM" button to change the current user.

Finally, press the "START/STOP" button to confirm user and turn off.


You can stop the setting anytime when you press the "START/STOP" button to save the current setting and turn the unit off.

If select to display the date and time, the current date and time will display when the unit is turned off.

17

Proper Measurement Method

- 1 Remove all clothes from your arm, allowing the cuff to fit directly on the skin.
- 2 Wrap the cuff tightly around your wrist by using the hook and loop strip.

- 3 Sit in a chair with your feet flat on the floor.


18

Take a Measurement

- 1 Press the "START/STOP" button.

All display symbols appear. If the positioning system is turned on, and the unit detects your arm is not in the proper position, it will sound the alarm to remind you to raise or lower the cuff. The arrowhead will then flash, then the arrowhead symbol will disappear, and the "♥" symbol will display.


19

Take a Measurement

- 2 When your arm is in the proper position and the heartbeat symbol disappears, the unit will automatically begin to inflate and measure.


The heartbeat symbol "♥" flashes once a pulse is detected.

- 3 The blood pressure and pulse rate are displayed when the measurement is complete.

The cuff deflates automatically and the measurement results are stored in memory with date and time.

The "♥" symbol will display if irregular heartbeat is detected.


20

Take a Measurement

4 Press the "START/STOP" button to turn the unit off.

The unit will automatically turn off after two minutes if you have forgotten to turn it off.

If selected to display the date and time, the current date and time will display when the unit is turned off.


NOTE:

To stop the inflation or measurement, press the "START/STOP" button at any time. The unit will stop inflating and then release and turn off.

21

Use the Memory Function

- The unit stores the blood pressure and pulse rate in the memory each time a measurement is completed. It can automatically store 180 sets of measurement values for 2 users (2 * 90). The earliest record is deleted to save the latest measurement value when any one user has more than 90 sets.
- The unit also calculates an average reading based on the values of the latest 3 measurements taken.
- ▶ To read the average value
- Press the "MEM" button; the average value of the current user's latest 3 measurements will display.


22

Use the Memory Function

▶ To read the measurement value

- Press the "MEM" button to display the set of measurement values with the measuring date from the most recent to the oldest.
- Press the "SET" button to display the next set of measurement values with the measuring date from the oldest to the most recent.


23

Use the Memory Function

▶ To delete the memory value

- In the memory reading mode, press the "SET" button for 4-6 seconds; the unit will display "EL".
- Press the "MEM" button and confirm "YES" or "NO".
- If you are sure all of the corresponding user's data can be deleted, press the "START/STOP" button; all the corresponding user's data will then be deleted and the unit will turn off.

Note: If you have set "YES", but do not press a button for 2 minutes, the unit will turn off automatically, but the corresponding user's blood pressure data will not be deleted. The delete operation must be confirmed by the user. After the measurement has finished and the memory values display, the backlight remains on for about 10 seconds.

24

About Blood Pressure

Blood Circulation


The blood circulation is responsible for supplying the body with oxygen.

Blood pressure is the pressure exerted on the arteries.


The systolic blood pressure value (higher pressure or top value) represents the blood pressure produced by contraction of the heart muscle.

The diastolic blood pressure value (lower pressure or lower value) represents the blood pressure produced by relaxation of the heart muscle.

Systolic Blood Pressure

Diastolic Blood Pressure


25

About Blood Pressure

Health and Blood Pressure

The incidence of hypertension increases with age. In addition, lack of exercise; excess body fat; and high levels of cholesterol (LDL), which sticks to the inside of blood vessels, reduce their elasticity. Hypertension accelerates arterial sclerosis, which can lead to very serious conditions such as stroke and myocardial infarction. For these reasons it is very important to know whether the blood pressure is within a healthy range. Blood pressure fluctuates from minute to minute throughout the day, making it essential to take regular measurements to help you identify an average blood pressure.


26

About Blood Pressure

Classification of Blood Pressure

After each measurement is completed, the LCD will automatically show your position on the six segments of the bar indicator which correspond to the WHO (World Health Organization) Blood Pressure Indicator, shown below.

Reference Material: Journal of Hypertension 1999, Vol 17 No.2

***Note!**

When a person's systolic and diastolic pressures fall into different categories, the higher category should apply.

27

About Blood Pressure

Classification of Blood Pressure

Symptoms of High Blood Pressure

High blood pressure can go unnoticed for a long time, since it doesn't cause noticeable symptoms. The following are all possible causes of abnormally high blood pressure:

- Overweight;
- High cholesterol level;
- Smoking;
- Excessive alcohol consumption;
- Stress and emotional upset;
- Excessive consumption of salt;
- Lack of physical exercise;
- Genetic / hereditary predisposition;
- Underlying illness, such as kidney disorders or endocrine disturbance.

Treatment of High Blood Pressure

If your blood pressure reaches upper values of 140~160 mmHg and lower values of 90~95 mmHg in repeated measurements over several days, you should consult your doctor for detailed medical examination. You can assist the treatment prescribed by your doctor in the following ways:

- Lose weight and lower your cholesterol level;
- Reduce the consumption of alcohol;
- Reduce the intake of salt;
- Stop smoking;
- Take regular exercise;
- Monitor your blood pressure.

28

Troubleshooting

Error Indicators

■ The following symbol will appear on the display when measuring abnormally.

Symbol	Cause	Correction
	Inflation error	Wrap the cuff correctly Inflate again, ensuring cuff is wrapped correctly
	Measurement has failed	Do not move your arm or body; remain still and quiet Measure again, ensuring correct procedure
	Battery power is too low	Replace both worn batteries with new ones

29

Troubleshooting

Eliminate Fault


Problem	Causes and Solutions
No power	Replace both worn batteries with new ones
Nothing appears on the display	Check the battery installation for proper placement of both battery polarities
Measurement values appear too high or too low	Blood pressure varies constantly. Many factors including stress, time of day, how you wrap the cuff, etc., may affect your blood pressure. Review the sections "Proper Measurement Method" and "Take a Measurement".

30

Care and Maintenance

Care

- Keep the unit in the storage case when not in use.
Clean the unit with a soft dry cloth.
- ▲ NOTICE: DO NOT use any abrasive or volatile cleaners.
- ▲ NOTICE: DO NOT immerse the unit or any components in water.

- ▲ NOTICE: Avoid subjecting the monitor to strong shocks or vibrations, such as dropping the unit on the floor.
- ▲ NOTICE: Use the unit in accordance with the instructions provided in this manual. Use only authorized parts and accessories.

31

Care and Maintenance

Maintenance

▲ NOTICE: DO NOT clean the body or cuff with solvents, naphtha, thinner or gasoline etc.


▲ NOTICE: DO NOT wet the cuff or attempt to clean the cuff with water.

▲ NOTICE: Store the unit in a clean, dry location. DO NOT subject the unit to extreme hot or cold temperature, humidity, or direct sunlight.

▲ NOTICE: Remove the batteries if the unit will not be used for 3 months or longer. Always replace both batteries with new ones at the same time.

32

Specifications

Model No.		Lumiscope 1146
Display		LCD Digital Display
Measuring Principle		Oscillometric Method
Measurable Circumference of Wrist		13.5~21.5cm (5.3"~8.5")
Electric Shock Protection		Internal power supply appliance type B
Backlight		White
Measurement Range	Pressure	0 to 299 mmHg (0~39.9 kPa)
	Pulse	40 to 180 Pulses/min
Accuracy	Pressure	± 3 mmHg (0.4kPa)
	Pulse	± 5% of reading
LCD Indication	Pressure	3 Digits Display of mmHg or kPa
	Pulse	3 Digits Display
	Symbol	Memory/IHB/Heartbeat/Low Battery, etc.

33

Specifications

Inflation		Automatic by Internal Pump
Rapid Air Release		Automatic by Air Valve
Memory Function		Memory 180 (2*90) Sets of Measurement Values
Power Source		2 AAA Alkaline Batteries
Automatic Power-Off		In 2 Minutes
Gross Weight		Approx. 170g (Batteries not included)
Operating Environment	Temperature	10 ~ 40 °C (50 ~ 104 °F)
	Humidity	15 ~ 90%RH (Noncondensing)
Storage Environment	Barometric Pressure	80~105 kPa
	Temperature	-20 ~ 60 °C (-4 ~ 140 °F)
	Humidity	10 ~ 95% RH (Noncondensing)

34

Specifications

- The unit satisfies the requirements of IEC60601-1 Medical electrical equipment, EN1060-1: Non-invasive sphygmomanometer, EN1060-3: Non-invasive sphygmomanometer, IEC80601-2-30: Medical electrical equipment.
- The unit satisfies the requirements of EN60601-1: Medical electrical equipment Part 1: General requirements for basic safety and essential performance; EN60601-1-2: Medical electrical equipment Part 1-2: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests; EN1060-1: Non-invasive sphygmomanometers – Part 1: General requirements; EN1060-3: Non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring system. EN1060-4: Non-invasive sphygmomanometers – Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers.

35

Specifications

- The unit has passed the EMC Test.
- Blood pressure measurements determined with this device are equivalent to those obtained by a trained observer using the cuff/stethoscope auscultatory method, within the limits prescribed by the American National Standard, Manual, electronic, or automated sphygmomanometers.
- The risk of patient and user can be lowered to an acceptable level.
- Specifications are subject to change without notice due to improvements in performance.
- Applied parts include cuff.
- The device, accessories, and packaging must be disposed of correctly at the end of their usage. Please follow Local Ordinances and Regulations for disposal.
- Protection Class: Internally powered equipment
- Applied Part Type: Type B
- Moisture Protection: IPX0, continue operation

36

Limited Warranty

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser ("Customer") only, that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use for the following period after the purchase date unless there is an expiration date on the component in which case the warranty shall expire on the earlier of the warranty period or the expiration date.

Blood Pressure Monitor: two years
Cuff: one year

The warranty does not extend to non-durable parts and does not include labor or costs of shipping. This limited warranty is not transferable. All warranties are conditioned upon the proper use of the product strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. The warranty is void if the defect is caused by any other reason not related to defects in materials or workmanship.

OBTAINING WARRANTY SERVICE

GF's customer service team must be notified of any warranty claim within the applicable warranty period. Call 770-368-4700, or fax 770-368-2386 or e-mail cs@grahamfield.com. There is a \$15 fee for return shipping and handling. Please make checks payable to GF Health Products, Inc. Failure to follow the specific directions provided by the GF customer service team will result in denial of the warranty claim.

37

Limited Warranty

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS.

GF'S TOTAL LIABILITY FOR ANY PRODUCT OR SERVICE PROVIDED IS LIMITED TO THE COST OF THE PRODUCT GIVING RISE TO THE CLAIM. IN NO EVENT WHETHER IN CONTRACT, INDEMNITY, WARRANTY, NEGLIGENCE, STRICT LIABILITY OR OTHERWISE WILL GF BE LIABLE FOR ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

This warranty gives you specific legal rights. You may have additional rights which vary from state to state.

38

Appendix A: EMC Information

- ▲ **NOTICE:** The use of ACCESSORIES, transducers and cables other than those specified, with the exception of transducers and cables sold by the MANUFACTURER of the Digital Blood Pressure Monitor as replacement parts for internal components, may result in increased EMISSIONS or decreased IMMUNITY of the ME EQUIPMENT or ME SYSTEM.
- ▲ **NOTICE:** This machine should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.
- ▲ **NOTICE:** MEDICAL ELECTRICAL EQUIPMENT needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS.
- ▲ **NOTICE:** Portable and mobile RF communications equipment can affect MEDICAL ELECTRICAL EQUIPMENT.
- ▲ **NOTICE:** Do not use a mobile phone, or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- ▲ **NOTICE:** This unit has been thoroughly tested and inspected to assure proper performance and operation.

39

Appendix A: EMC Information

Guidance and manufacturer's declaration

Guidance and manufacturer's declaration – electromagnetic emission		
The Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer of the user of the Digital Blood Pressure Monitor should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Digital Blood Pressure Monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The Digital Blood Pressure Monitor is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

40

Appendix A: EMC Information

Guidance and manufacturer's declaration – electromagnetic immunity			
The Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Digital Blood Pressure Monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floor are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

41

Appendix A: EMC Information

Guidance and manufacturer's declaration – electromagnetic immunity			
The Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Digital Blood Pressure Monitor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3A/m 80 MHz to 2.5 GHz	3A/m	Portable and mobile RF communications equipment should be used no closer to any part of the Digital Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.167\sqrt{P}$ 80 MHz to 800 MHz $d = 2.333\sqrt{P}$ 800 MHz to 2.5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site.

42

Appendix A: EMC Information

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Radiated RF IEC 61000-4-3	3A/m 80 MHz to 2.5 GHz	3A/m	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Digital Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Digital Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Digital Blood Pressure Monitor.			

43

Appendix A: EMC Information

Recommended separation distances between portable and mobile RF communications equipment and the Digital Blood Pressure Monitor.		
The Digital Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Digital Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Digital Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment.		
Output power of transmitter (W)	Separation distance according to frequency of transmitter (m)	
	80 MHz to 800 MHz $d = 1.167\sqrt{P}$	800 MHz to 800 MHz $d = 2.333\sqrt{P}$
0.01	0.117	0.233
0.1	0.369	0.738
1	1.167	2.333
10	3.689	7.379
100	11.667	23.333
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.		

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Graham-Field, Lumiscope, y Lumiscope For The Quality Of Life son marcas registradas de GF Health Products, Inc.
 © 2014 GF Health Products, Inc. Fabricado en China

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