JOHN BUNN®
by graham-field

DigiO₂

Finger Pulse Oximeter

Operation Manual

Important: Do not operate the John Bunn DigiO₂ Finger Pulse Oximeter without first reading and understanding this manual! Save this manual for future use.

JBO2007-INS-LAB-RevD18

CONTENTS

INTRODUCTION	3
IMPORTANT SAFETY GUIDELINES - PLEASE READ BEFORE USE	3
INTENDED USE	6
CONTRAINDICATIONS	6
DECLARATION OF CONFORMITY:	6
OVERVIEW	7
APPEARANCE	7
FEATURES	8
NAME AND MODEL	8
BATTERY INSTALLATION	8
OPERATION INSTRUCTIONS	9
ACCESSORIES (INCLUDED)	12
REPAIR AND MAINTENANCE	12
MAINTENANCE	12
CLEANING AND DISINFECTION	13
TROUBLESHOOTING	14
SYMBOL KEY	15
TECHNICAL SPECIFICATIONS	15
LIMITED WARRANTY	17
DIGIO ₂ OXÍMETRO DE PULSO MANUAL DE OPERACIÓN (EN ESPAÑOL)	19

INTRODUCTION

Dear User: Please read and understand all instructions before using the John Bunn DigiO₂ Finger Pulse Oximeter, JBO2007.

Important Safety Guidelines - Please Read Before Use

The safety statements presented in this chapter refer to the basic safety information that the operator of the John Bunn DigiO₂ Finger Pulse Oximeter shall pay attention to and abide by. There are additional safety statements in other chapters or sections, which may be the same as or similar to the following, or specific to the operations.

- ⚠ WARNING:** Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious injury.
- ⚠ CAUTION:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor or moderate personal injury.
- ▲ NOTICE:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor personal injury or product/property damage.

Info: Provides application recommendations or other useful information to ensure that you get the most from your product.

- ⚠ WARNING: Important! Read and understand this manual before using the oximeter. Do not use this device without proper instruction from a Healthcare Professional. If the oximeter is not properly operated, personal injury and/or damage to the oximeter could result.**

- ⚠ WARNING: If components are damaged or missing, contact your dealer immediately. DO NOT use substitute parts.**
- ⚠ WARNING: Inspect the device before use to ensure that there is no visible damage that may affect patient or user's safety or measurement performance. Thereafter, inspect the device at least once a week.**
- ⚠ WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper installation or use of this product.**
- ⚠ WARNING: DO NOT open the oximeter housing. It contains no serviceable parts. Service must be performed by trained, authorized personnel only. Otherwise, device failure and health hazard may occur.**
- ⚠ WARNING: EXPLOSION HAZARD: Do not use the oximeter in the presence of flammable anesthetics, explosive substances, vapors or liquids.**
- ⚠ WARNING: The oximeter is designed for real-time and rapid measurement. It is not suitable for long-time continuous patient monitoring. DO NOT apply the sensor to the same finger for more than two hours. If any readings are found to be abnormal, please change the position of the sensor.**
- ⚠ WARNING: The oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with clinical signs and symptoms. It is intended as a measuring, not a treatment, device.**

- ⚠ WARNING: Magnetic and electrical fields are capable of interfering with the proper performance of the oximeter. For this reason, ensure that all external devices operated in the vicinity of the oximeter comply with the relevant EMC requirements. Mobile phones, X-ray equipment, and/or MRI devices are a possible source of interference, as they may emit high levels of electromagnetic radiation.**
- ⚠ WARNING: DO NOT clip this device on edema or tender tissue.**
- ⚠ WARNING: The infrared light emitted from the device (it is invisible) is harmful to the eyes. DO NOT stare into the light.**
- ⚠ WARNING: Keep this device out of the reach of children.**
- ⚠ WARNING: Cancer and Reproductive Harm - www.p65warnings.ca.gov.**
- ⚠ WARNING: To ensure patient safety, use only parts and accessories specified in this manual.**
- ⚠ WARNING: Keep the oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.**
- ⚠ WARNING: If the oximeter gets wet, stop operation immediately and do not resume operation until it is dry. When it is carried from a cold environment to a warm and/or humid environment, allow it to warm to room temperature before using it.**
- ▲ NOTICE: DO NOT use sharp objects to press the front panel button.**
- ▲ NOTICE: DO NOT immerse the oximeter in liquid.**

- ▲ NOTICE: DO NOT** subject the oximeter to high temperature or high pressure steam disinfection. Refer to the **REPAIR AND MAINTENANCE/Cleaning and Disinfection** section of this manual for instructions for cleaning and disinfection.

Intended use

The John Bunn DigiO₂ Finger Pulse Oximeter, which combines a sensor and monitor into one device, is used to measure pulse oxygen saturation (SpO₂) value, pulse rate value, and blood perfusion index (PI) or pulse strength in home or clinical settings. It is designed for adult and pediatric patients, and fits fingers of thickness .39" - .87" (10mm - 22mm).

Contraindications

- ⚠ WARNING: DO NOT** use the oximeter while the patient is undergoing MRI or CT scanning.

Declaration of Conformity:

This device complies with the following standards:	IEC 60601-1
	ISO 9919
	ISO 21647
	and follows the provisions of the council directive MDD93/42/EEC

OVERVIEW

Appearance


Figure 1
front view


Figure 2
rear view

Features

The John Bunn DigiO₂ Finger Pulse Oximeter provides a simple way to spot-check patients by combining the sensor and monitor into one integrated, compact, easy to use device with single-button operation. The oximeter can accurately measure pulse oxygen saturation (SpO₂) value, pulse rate value, and blood perfusion index (PI) or pulse strength. When a finger is inserted into the sensor's rubber cushion, the SpO₂ value and pulse rate value automatically display on the screen. The oximeter features a bright-color LED screen with multi-directional displays that enable the user to read results from any angle. The oximeter features both audible and visible alarms and a low voltage indicator. The oximeter powers off automatically in eight seconds when not in use. The oximeter provides approximately 2,000 spot checks with the use of alkaline batteries.

Name and Model

Name: John Bunn DigiO₂ Finger Pulse Oximeter

Model: JBO2007

BATTERY INSTALLATION

1. Refer to Figure 3; insert two AAA batteries, oriented as shown, into the battery compartment.
2. Replace the cover.

Info: Ensure the batteries are correctly installed. Incorrect installation will prevent the device from operating.


Figure 3
battery installation

OPERATION INSTRUCTIONS

⚠ CAUTION: Do not use the sensor if it is damaged.

Info: Ensure the finger is dry before being inserted into the oximeter.

Info: Do not shake the finger; keep finger motionless during oximeter use.

Info: Avoid placing the device on a limb which is wrapped with a cuff for blood pressure measurement or undergoing venous infusion.

Info: Ensure that no obstruction (contamination, scar tissue, fingernail polish, acrylic fingernail, etc.) exists at the site where the sensor is placed; such obstruction could affect the signal received by the sensor and result in an incorrect measurement.

Info: Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.

1. Open the clip as shown in Figure 4.
2. Insert finger into the clip's rubber cushion (ensure that the finger is correctly positioned), and then gently clip the finger.
3. The device will power on automatically in two seconds, and start to display the software version number.

**Figure 4
operation**

4. The screen will then display the data as shown in Figure 5. After the waveform and value become stable, the user can read the information on the display screen.
5. When the screen displays as shown in Figure 5, press the Display Key to change the display screen's orientation:
 - a. Press the Display Key once; the display screen will rotate 180° and display as shown in Figure 6.
 - b. Press the Display Key twice; the screen will now display as shown in Figure 7.
 - c. Press the Display Key three times; the display screen will rotate 180° and display as shown in Figure 8.
 - d. Press the Display Key four times; the display screen will return to the original screen as shown in Figure 5.

Figure 5

Figure 6

Figure 7

Figure 8

6. Press the Display Key for approximately two seconds; the screen will display as shown in Figure 9. Differences between Figure 9 and Figure 5 are as follows:
 - a. In Figure 5, **SpO₂** and **PR (pulse rate)** are being monitored and displayed on the screen.
 - b. In Figure 9, **SpO₂** and **PI (blood perfusion index)** are being monitored and displayed on the screen.

Figure 9

Figure 10

7. Press the Display Key twice; the screen will now display as shown in Figure 10.
8. The Display screen as shown in Figure 9 or Figure 10 will return to the original screen as shown in Figure 5 or Figure 6 if without operation in ten seconds.
9. Alarm Indicator: When measuring, if SpO₂ value and pulse rate value are outside the preset alarm limits, the device will alarm automatically and the value which is outside the limit will flash on the screen; press the Display Key to suspend the alarm beep (see alarm table below).

Alarm	SpO ₂	Outside limit alarm: three beeps	Preset limit: ≤90%
	Pulse Rate	Outside limit alarm: two beeps	Preset limits: Low ≤50bpm, High ≥120bpm

ACCESSORIES (INCLUDED)

1	One lanyard
2	Two AAA alkaline batteries
3	One pouch
4	One operation manual
5	One quality certificate

REPAIR AND MAINTENANCE

⚠ WARNING: DO NOT open the oximeter housing. It contains no serviceable parts. Service must be performed by trained, authorized personnel only. Otherwise, device failure and health hazard may occur.

Maintenance

In order to ensure the oximeter's longevity and continued dependable operation, please observe the following maintenance recommendations.

Change the batteries when the low-voltage indicator displays.

Clean the surface of the device before using: Wipe the device with alcohol, then let it air dry or wipe it dry. Repeat as necessary.

Take out the batteries if the oximeter will not be used for a long period of time.

Keep the device in a dry environment. Humidity may shorten the life of the device, or even damage it.

Environmental requirements	Operation	Temperature: 41°F ~104°F (5°C ~40°C)
		Relative humidity: 30%~80%
		Atmospheric pressure: 10.15 psi~15.37 psi (70kPa~106kPa)
	Storage	Temperature: -4°F to 131°F (-20°C to 55°C)
Relative humidity: < 95%		

Cleaning and Disinfection

1. Do not immerse the device in liquid. Clean the device only by using a surface-clean technique: Wipe with a clean, soft cloth saturated with a solution of 75% isopropyl alcohol.

If low-level disinfection is required, wipe with a 1:10 bleach to water solution.

2. Wipe with a clean, soft cloth saturated with clean water and air dry, or dry with a clean, soft cloth.

▲ **NOTICE: Do not immerse the device in liquid.**

▲ **NOTICE: Do not clean the device with abrasive cleansers.**

▲ **NOTICE: Do not sterilize the device with irradiation, steam, or ethylene oxide.**

▲ **NOTICE: Do not disinfect the device with high temperature or high pressure steam.**

TROUBLESHOOTING

Problem	Possible Reason	Solution
The SpO ₂ and Pulse Rate display is unstable	1. The finger is not correctly positioned in the device	1. Position the finger correctly and try again
	2. The finger is shaking or the patient is moving	2. Ensure that patient is able to remain immobile
The device does not turn on	1. The batteries are low	1. Replace the batteries
	2. The batteries are not correctly installed	2. Reinstall batteries
	3. The device has malfunctioned	3. Contact your GF distributor
The indicator light suddenly turns off	1. The device powers off automatically when it gets no signal for eight seconds	1. Normal
	2. The batteries are almost drained	2. Replace the batteries

SYMBOL KEY

Symbol	Description
	Type BF
	WARNING: See operation manual
%SpO₂	The pulse oxygen saturation
PI	Perfusion Index or Pulse Strength
 bpm	Pulse rate (beats per minute)
	Low battery voltage
SN	Serial number
	Do not treat this device as household waste—separate collection required

TECHNICAL SPECIFICATIONS

Display mode	LED Display
Power supply requirement	Two 1.5V AAA alkaline batteries
	Supply voltage: 2.7V~3.3V
Operating current	≤40mA

TECHNICAL SPECIFICATIONS continued	
SpO ₂ parameter specifications	Measuring range: 35%~99%
	Accuracy: ±2% during 75%~99%, ±3% during 50%~75%
Pulse rate parameter specifications	Measuring range: 25-250 bpm
	Accuracy: ±2 bpm or ±2%
Pulse strength parameter specifications:	Measuring range: 0.2-20%
	Accuracy: ±0.1% during 0.2%~2%, ±1% during 2%~10%, ±2% during 10%~20%
Resistance to interference of surrounding light	The difference between the value measured in the condition of indoor natural light and that of darkroom is less than ±1%
Resistance to interference of man-made light	Values of SpO ₂ and pulse rate can be accurately measured by pulse oxygen simulator
Dimensions (L x W x H)	2.59" x 1.42" x 1.3" (66 mm x 36 mm x 33 mm)
Net weight with batteries	2 oz. (60g)
Classification	The type of protection against electric shock: Internally powered equipment
	The degree of protection against electric shock: Type BF applied part
	The degree of protection against harmful ingress of liquids: Ordinary equipment without protection against ingress of water
	Electromagnetic compatibility: Group I, Class B

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted. This limited warranty shall only apply to defects that are reported to GF's customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable. Within the guidelines set forth in this document, this product is warranted for three (3) years. The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

This limited warranty shall only apply to defects that are reported to the Distributor from whom the Customer purchased the product within the applicable warranty period. If there is not a Distributor, you must contact GF directly by calling 1.770.368.4700, sending a fax request to 1.770.368.2386, or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: tips, casters, filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 6) Any labor or shipping charges incurred in the replacement part installation or repair;
- 7) Costs and expenses of regular maintenance and cleaning; and
- 8) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.

1.770.368.4700

Information contained herein is subject to change.

The most current and complete product information can be found on our website.

www.grahamfield.com

Follow us on

© 2007, GF Health Products, Inc. All Rights Reserved.

John Bunn and Graham-Field are trademarks of GF Health Products, Inc.,

One Graham-Field Way, Atlanta GA 30340-3140

GF Health Products, Inc. is an ISO 13485 Certified Company.

Manufactured for GF Health Products, Inc.

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John Bunn y Graham-Field son marcas de GF Health Products, Inc.
GF Health Products, Inc. ha sido certificado por la norma ISO 13485.
Fabricado para GF Health Products, Inc.