
LUMISCOPE®

by graham-field

AUTOMATIC DIGITAL BLOOD PRESSURE MONITOR

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
SYMBOLS

Symbols	Meaning
	<p>Symbol for the marking of electrical and electronics devices according to Directive 2002/96/EC.</p> <p>The device, accessories and the packaging have to be disposed of waste correctly at the end of the usage. Please follow Local Ordinances or Regulations for disposal.</p>
	<p>Keep dry</p>
	<p>Attention, consult accompanying documents</p>
	<p>Type BF Applied Part</p>
	<p>standby</p>

TABLE OF CONTENTS

1	GENERAL	5
	PRINCIPLE OF OPERATION	5
	NEW TECHNOLOGIES USED	5
2	PREPARATION FOR MEASUREMENT	6
3	BATTERY INSTALLATION	8
4	CORRECT POSTURE FOR MEASUREMENT	8
5	ASSEMBLING THE CUFF	8
6	SETTING THE DATE AND TIME	10
7	CARRYING OUT A MEASUREMENT	11
	IRREGULAR HEARTBEAT DETECTION	12
8	MEMORY FUNCTION	13
	MEMORY RECALL	13
	MEMORY CLEARANCE	13
9	WHO BLOOD PRESSURE CLASSIFICATION INDICATION	14
11	ERROR AND LOW BATTERY INFORMATION	15
12	MAINTENANCE, STORE, REPAIR AND RECYCLE	15
13	TROUBLESHOOTING	16
14	SPECIFICATIONS	17
15	LIMITED WARRANTY	18
16	BLOOD PRESSURE RECORD	19
17	MANUFACTURER'S DECLARATION	20

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GENERAL

Intended Use: This device is intended for the non-invasive measurement of systolic and diastolic arterial blood pressure and pulse rate in adults (age 15 and above).

Contraindication: Consult a physician if measurement is taken in children or persons with arrhythmia, as errors may occur.

This instruction manual is intended to assist the user for safe and efficient operation of the automatic digital talking blood pressure monitor (hereinafter: device) model 1133. The device must be used in accordance with the procedures described in the manual. Read and understand the entire manual, especially the section **<Tips for taking blood pressure measurements>**.

PRINCIPLE OF OPERATION

This device adopts the oscillometric technology with Fuzzy Algorithm to measure the arterial blood pressure and pulse rate. The cuff is wrapped around the arm and automatically inflated by the air pump. The sensor of the device catches weak fluctuation of the pressure in the cuff produced by extension and contraction of the artery of the arm in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters of the mercury column, and is displayed by digital value.

⚠ WARNING: This device cannot provide reasonable accuracy if used or stored in the temperature or humidity beyond the range stated in the <SPECIFICATIONS> section of this manual.

⚠ WARNING: DO NOT USE THE DEVICE OUTDOORS

TECHNOLOGIES USED

Fuzzy Algorithm is the processing algorithm, taking into account the specialty of individual heartbeats, which provides higher accuracy of measurement.

Proposition 65 Warning

⚠ WARNING: Cancer and Reproductive Harm - www.p65warnings.ca.gov.

PREPARATION FOR MEASUREMENT

1. Do not conduct any measurements if the temperature is low (below 50°F / 10°C) or high (over 104°F / 40°C), or if the relative humidity is higher than 85%, as this can lead to inaccurate readings.
2. Take the measurement at room temperature in a quiet and stress-free environment.
3. Do not move yourself or the unit during the measurement. Do not speak during the measurement.
4. Blood pressure varies naturally depending on the time of day and is affected by many factors. Blood pressure is usually highest at work and reaches its lowest level during sleep.
5. Blood pressure measurements should be assessed by a physician or trained healthcare professional who is familiar with your medical history. If you use the unit and regularly record the results, keep your physician informed of any changes in your blood pressure.
6. The performance of this device can be affected as severe arrhythmias such as atrial or ventricular premature beats or atrial fibrillation are presented during measurement.
7. The blood pressure measurements conducted with this unit are equivalent to measurements obtained by a trained observer in accordance with the values achieved using the cuff/stethoscope auscultation method and are within the specified EN 1060-4 standard limits.
8. Patients with cardiovascular disease: take measurement under physician's instruction. Under no circumstances should you alter the dosage of any drug prescribed by your physician!
9. Accurate measurement of blood pressure may be difficult in case of arrhythmia, premature beat, atrial fibrillation, arteriosclerosis, hypoperfusion, diabetes, pregnancy, nephropathy, weak pulse, or in patients with obvious fluctuation of heart contraction rhythm. Consult a physician to interpret your blood pressure readings.

10. It is necessary to keep quiet during the measurement to get an accurate reading. Measurement should be conducted in a quiet environment at room temperature. Don't eat or smoke before a measurement.

This device is supplied with the standard cuff which fits arm sizes 8.7-13.8 in. / 22-35 cm.

⚠ WARNING: Ensure the cuff size is appropriate for the person whose blood pressure is being taken.

⚠ WARNING: Do not use any cuff other than the original cuff contained in this kit!

11. Repeated measurements with an interval of 3 minutes are recommended, so you can calculate the average to get a more accurate measurement. Atherosclerosis patients may require a longer interval (10-15 minutes) as elasticity of the patient's vessels decreases significantly with the disease. 10-15 minute interval is also applicable for patients suffering from longtime diabetes.

CLASSICATION

- ME EQUIPMENT is not intended for use in an oxygen rich environment or in the presence of flammable mixers.
- Internally powered equipment.
- Type BF applied part, cuff recognized as applied part.

BATTERY INSTALLATION

1. Open the battery cover and install four AA alkaline batteries in the battery compartment as indicated. Ensure the polarity is correct;
2. Close the battery compartment cover.
 - Replace the batteries when the replacement indication “” appears in the display or nothing after “” button is pressed;
 - Do not use rechargeable batteries;
 - Replace all batteries simultaneously;
 - If the device will be unused for a long period of time, remove the batteries;
 - DO NOT leave worn batteries in the device.

CORRECT POSTURE FOR MEASUREMENT

1. Sit beside the table and let the table support your arm as you take the measurement. Ensure the cuff on the upper arm is at approximately the same level as the heart as shown in Fig.1, and ensure your feet lie flat on the ground and do not cross.


Fig.1

ASSEMBLING THE CUFF

1. Insert the edge of the cuff approximately 2 in. / 5 cm into the D-ring as shown in Fig. 2.







Fig.2

- Put the cuff on the left upper arm with the tube pointing toward the palm. If measurement on your left arm is difficult, you may use the right arm for measurement. In this case, be aware that the readings may differ about 5-10 mmHg between left arm and right arm.
- Wrap cuff around your upper arm with the lower edge of the cuff approximately .8 - 1.2 in. / 2-3 cm above the elbow. The mark <ARTERY> must be over the artery of the upper arm.
- Press the cuff to ensure it is attached securely. The cuff should not be too tight or loose. Two fingers should easily fit between cuff and upper arm.
- The mark <INDEX> on the cuff must point to area <NORMAL> or <LARGE CUFF>. This means the cuff size is correct. If mark <INDEX> points to the area beyond area <NORMAL> or <LARGE CUFF>, consult your distributor to see whether you need another size cuff.
- Sometimes it is difficult to make the cuff fit well, depending on the shape of the user's upper arm; the cone-shape assembly of cuff is also acceptable.
- If your clothes restrict the blood circulation of your upper arm, or you roll your sleeve up so as to result in such restriction, take off your shirt to get an accurate measurement if necessary.

SETTING THE DATE AND TIME

This function provides accurate measuring time for each measurement. To get an accurate date and time, preset the date and time correctly before the first use of this device.

The operation procedure for presetting Date/Time is as follows:

1. When the device is connected to power supply at first, the display will show as in Fig.3;
- 2 Press and hold the button 'M1', then press the button '⏻', and the year number flashes as shown in Fig.4;
3. Press button 'M1', or 'M2' to subtract or add the number, and press button '⏻' for confirmation;
4. When the year setup is finished, the month number will flash automatically. Follow the same instruction as above to set month, date and time;
5. Press button '⏻' to finish the setup. If you wish to change the date and time, repeat steps 2,3,4.

Fig.3



Fig.4

Note: When under the mode of Date/Time setting without any operation, the device will automatically return to standby mode after one minute.


CARRYING OUT A MEASUREMENT

1. Assemble the cuff, then insert the tube plug into the air connector. Before the measurement, take 3~5 deep breaths and relax. Don't talk or move your arm;
2. Press button '⏸', and all symbols will appear on display in 2 seconds then '0 mmHg' will appear on the screen. The pump begins to inflate with display showing the pressure reading. Generally the pressure will reach 190 mmHg as shown in Fig.5;


3. The pump stops inflating and the pressure begins to decrease gradually, during which the user's blood pressure and pulse will be calculated as shown in Fig.6;

IRREGULAR HEARTBEAT DETECTION

When the device detects the irregular heartbeat or any excessive body movement during measurement, '♥' icon will flash on the screen as shown in Fig.7.

⚠ WARNING: Contact your physician if you see the '♥' indicator frequently.

- The cuff will deflate quickly and the blood pressure reading, pulse reading will show in the display. Then the icon '👤' will flash to remind the user to record the reading as shown in Fig 8;

Fig.8

- Press button 'M1' or button 'M2' to record the reading in the corresponding memory; the measurement time will be displayed on the screen alternately. If the user does not press the button, the reading won't be recorded;
- Press the button '⏻' to turn off the device. Rest for at least 5 minutes before taking another measurement. If the power supply is not switched off and the device is unused for 3 minutes, the device will switch off automatically.

AUTOMATIC INFLATION

There are four inflation pressure levels for this device: 190mmHg, 230mmHg, 270mmHg and 300mmHg. When 190mmHg is not enough or arm movement occurs, the device will automatically inflate to higher pressure to ensure a successful measurement. It is not a fault.

RAPID DEFLATION DURING MEASUREMENT

If you do not feel well during measurement or want to stop the measurement for some reason, press the button '⏻'. The device will quickly release the air in the cuff and the device will be returned to standby mode.

MEMORY FUNCTION

MEMORY RECALL

1. The device can store 90 sets of readings respectively in 'M1' and 'M2', and will automatically calculate the average value of the latest 3 readings for 'M1' and 'M2'. When the memory is full (90 sets of readings are stored), the oldest reading will be replaced by a new one. Readings in memory will not clear even if power supply is removed;
2. When a measurement is finished or the device is on standby, the user can press button 'M1' or 'M2' to recall memory. Press button 'M1' or 'M2', the display will show the average value of the latest 3 readings as shown in Fig.9;
3. Press again, the display will show '01', which means the latest reading, then turns to another screen to show readings and measuring time as shown in Fig.10;
4. Press again, the display will show '02', which means the second to the last reading.

Fig.9

Fig.10

MEMORY CLEARANCE

When a measurement is finished or when the device is on standby, hold down button 'M1' or 'M2' for at least 5 seconds, the display will show 'CLR' which means the stored reading for 'M1' or 'M2' is removed as shown in Fig.11.


Fig.11

WHO BLOOD PRESSURE CLASSIFICATION INDICATION

Standards for assessment of high or low blood pressure, regardless of age, have been established by World Health Organization (WHO) as shown in the chart below:

The indicator displays a segment, based on the current data, corresponding to the WHO classification.

For example, if your blood pressure is 145mmHg (systolic pressure), 88mmHg (diastolic pressure), according to the World Health Organization standard, your blood pressure level is Mild Hypertension.

Note: If the systolic blood pressure and the diastolic blood pressure fall into different categories, the higher value should be taken for classification.

ERROR AND LOW BATTERY INFORMATION

INDICATION	POSSIBLE REASON	CORRECTION METHODS
	<p>The cuff is applied incorrectly or the tube plug is inserted loosely</p> <p>Arm/hand movement or talking occurred during measurement</p> <p>Arrhythmia</p>	<p>Ensure cuff is correctly applied and tube plug is inserted tightly and repeat the measurement</p> <p>Follow the recommendations in this manual precisely and repeat the measurement</p> <p>Consult your physician</p>
	<p>Weak batteries</p>	<p>Replace all four batteries with new ones</p>

MAINTENANCE, STORING, REPAIR AND RECYCLING

1. Protect this device against moisture, direct sunlight, shock, solvent, alcohol and gasoline.
 2. Remove the batteries if the device is being stored for a long period of time.
- ⚠ WARNING: Keep the device and batteries away from children.**
3. Keep sharp objects away from the cuff. Do not extend or twist the cuff.
 4. Use only soft, clean cloth to clean the device.
 5. The cuff is sensitive and must be handled with care. Clean the cuff with a clean, damp cloth for daily maintenance.
 6. To avoid cross-infection when sharing the cuff, sterilize the fabric cover of the cuff with a soft, clean cloth moistened with a 3% solution of hydrogen peroxide. After long use there will be a partial discoloration on the fabric surface of the cuff. Do not launder or iron the cuff.
- ▲ NOTICE: DO NOT wash the inner bladder.**
7. Follow your local recycling rules and dispose of device and batteries at an appropriate collection site.
 8. We recommend having the device inspected every 2 years to ensure proper function, accuracy and safety. Contact your distributor for maintenance.
- ⚠ WARNING: Do not open or repair the device. Contact your distributor for maintenance.**

TROUBLESHOOTING

SYMPTOM	CHECK POINT	REMEDY
No display after installing batteries	Depleted batteries Battery polarity incorrect Dirty battery compartment contact	Replace all batteries with new ones Correct battery polarity Clean battery terminals with dry cloth
Inflation stops and starts	Automatic inflation ensures correct measurement Talking or moving arm or hand during measurement	See<AUTOMATIC INFLATION>. Keep quiet and silent during measurement
The reading is extremely low or high	Is the cuff at the same level as the heart? Cuff wrapped incorrectly Moving arm during measurement Talking or moving arm or hand during measurement	Ensure correct posture Wrap the cuff correctly Relax during measurement Keep quiet and silent during measurement
Pulse rate is too low or too high	Did you talk or move your arm or hand during measurement? Did you take measurement right after exercise?	Keep quiet and silent during measurement Take measurement again after resting for at least 5 minutes
Batteries deplete quickly	Faulty batteries	Use new alkaline batteries of known manufacturers
Device automatically turns off	This is the result of automatic turn-off feature	This feature saves power consumption of the device, and it is not a fault

SPECIFICATIONS

Model	1133
Size (dimensions LxWxH)	4.13 x 4.72 x 2.17 in. / 105 ×120 ×55 mm
Weight (without batteries)	Approximately 11 oz. / 310 g
Measuring method	Oscillometry
Extreme Pressure	300mmHg
Measuring range	40 to 260 mmHg(blood pressure) 40 to 160 beats/minute(pulse rate)
Measuring accuracy	±3mmHg for static pressure ±5% of the reading for the pulse rate
Inflation	Automatic by pump
Rapid deflation	Automatic electronic valve
Batteries	Four AA 1.5V (not included)
Memory	2 Users with 90 sets of memory each
Operation temperature and humidity, air pressure	50°F to 104°F / 10°C to 40°C, ≤85% RH, 700hPa to 1060hPa
Transport and storage temperature and humidity, air pressure	-4°F to 122°F / -20°C to 50°C, ≤85% RH, 500hPa to 1060hPa
Upper arm circumference	fits arm circumference 8.7-13.8 in. / 22-35 cm (standard cuff)
Complete kit	Main body, cuff, gift box, instruction manual
Pollution Degrees	Degrees 2
Overvoltage category	Category II
Altitude	≤6,652 ft / 2,000 m

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. (See Obtaining Warranty Service below) This limited warranty is not transferable.

The warranted components and time periods are set forth below:

Digital Blood Pressure Monitor.....	24 months
Cuff.....	12 months

* Labor is not included in the warranty.

† Upholstery is only warranted on material supplied by GF.

‡ The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified.

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1.770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim.

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Matching of color, grain or texture except to commercially acceptable standards;
- 6) Changes in color caused by natural or artificial light;
- 7) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 8) Any labor or shipping charges incurred in the replacement part installation or repair;
- 9) Costs and expenses of regular maintenance and cleaning; and
- 10) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.

BLOOD PRESSURE RECORD

DATE	SYSTOLIC (mmHg)	DIASTOLIC (mmHg)	PULSE (beats/minute)

MANUFACTURER'S DECLARATION

The model 1133 Digital Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the model 1133 Digital Blood Pressure Monitor should assure that is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 6000-4-2	±6kV contact ±8kV air	±6 kV contact ±8 kV air	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in UT) for 0.5 cycle 40% U_T (60% dip in UT) for 5 cycles 70% U_T (30% dip in UT) for 25 cycles <5% U_T (>95% dip in UT) for 5 sec	<5% U_T (>95% dip in UT) for 0.5 cycle 40% U_T (60% dip in UT) for 5 cycles 70% U_T (30% dip in UT) for 25 cycles <5% U_T (>95% dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the 1133 Digital Blood Pressure Monitor Equipment requires continued operation during power mains interruptions, it is recommended that the 1133 Digital Blood Pressure Monitor Equipment be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3Vrms 150 kHz to 80 MHz	3 Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of the 1133 Digital Blood Pressure Monitor Equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = \left[\frac{3.5}{V_t} \right] \sqrt{P}$ $d = \left[\frac{3.5}{E_t} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{7}{E_t} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^c should be less than the compliance level in each frequency range.^d Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note 1 At 80MHz and 800MHz, the higher frequency range applies.

Note 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strength from fixed transmitters, such as base stations for radio telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 1133 Digital Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Model 1133 Digital Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model 1133 Digital Blood Pressure Monitor.

^b Over the frequency range 150KHz to 80MHz, field strength should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the 1133 Digital Blood Pressure Monitor

The 1133 Digital Blood Pressure Monitor is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the 1133 Digital Blood Pressure Monitor can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the 1133 Digital Blood Pressure Monitor as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{V_f} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_f} \right] \sqrt{P}$	$d = \left[\frac{7}{E_f} \right] \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

LUMISCOPE®

by graham-field

Manufactured For:
GF Health Products, Inc.

One Graham-Field Way, Atlanta, GA 30340

•1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

© 2008, GF Health Products, Inc. All Rights Reserved.
Lumiscopé is a registered trademark of GF Health Products, Inc.

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© 2008, GF Health Products, Inc. Todos los derechos reservados.
Lumiscope es una marca de GF Health Products, Inc.