WRIST DIGITAL BLOOD PRESSURE MONITOR

MODEL LUMISCOPE® 1143

INSTRUCTION MANUAL
1143-INS-LAB-RevF23

CONTENTS

INTRODUCTION..... 1
 IMPORTANT SAFETY PRECAUTIONS 2
 PRODUCT DESCRIPTION AND FEATURES 3
 SETUP – PREPARATION FOR OPERATION 4
 OPERATION..... 5
 MEASUREMENT INTERPRETATION 8
 TROUBLESHOOTING 9
 MAINTENANCE 10
 SPECIFICATIONS 11
 LIMITED WARRANTY 11
 MANUAL DE INSTRUCCIONES (ESPAÑOL)..... 13

INTRODUCTION

INTENDED USE

This device is intended for the noninvasive measurement of systolic and diastolic arterial blood pressure and pulse rate in adults. The device is intended to be used by adults, at room temperature, with skin contact. The patient is the intended operator.

Contraindications

- ⚠ **WARNING: DO NOT use this device on infants, children, or those who cannot express their own intentions.**
- ⚠ **WARNING: Consult your physician before use during pregnancy.**
- ⚠ **WARNING: Consult your physician before use if diagnosed with arrhythmia, arteriosclerosis, cardiovascular disease, diabetes, liver or kidney disease, severe hypertension, or peripheral circulatory disorders.**
- ⚠ **WARNING: This device is not AP/APG equipment and is not suitable for use in the presence of flammable gas (such as anesthetic gas, oxygen or hydrogen) or flammable liquid (such as alcohol).**

The device must be used in accordance with the procedures described in this manual. Read and understand the entire manual before use, especially the Important Safety Precautions section.

PRINCIPLE OF OPERATION

This device uses oscillometric technology with Fuzzy Algorithm to measure the arterial blood pressure and pulse rate. The wrist cuff is wrapped around the wrist and automatically inflated by an air pump. The sensor of the device catches weak fluctuation of the pressure in the wrist cuff produced by extension and contraction of the artery of the wrist in response to each heartbeat. The amplitude of the pressure waves is measured, converted in millimeters of the mercury column, and is displayed by digital value.

Info: This device cannot provide reasonable accuracy if used or stored at temperature or humidity beyond the range stated in the Specifications section of this manual.

IMPORTANT SAFETY PRECAUTIONS – READ BEFORE USE

The safety statements presented in this chapter refer to the basic safety information that the Blood Pressure Monitor (BPM) user shall pay attention to and abide by. There are additional safety statements in other chapters or sections, which may be the same as or similar to the following, or specific to the operations. Please note the following special statements, used throughout this manual, and their significance:

- ⚠ **WARNING: Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious personal injury.**
- ⚠ **CAUTION: Indicates a potential hazard or unsafe practice that, if not avoided, could result in minor or moderate personal injury.**
- ⚠ **NOTICE: Indicates a potential hazard or unsafe practice that, if not avoided, could result in product or property damage.**

Info: Provides application recommendations or other useful information to ensure that you get the most from your product.

⚠ **WARNING: Important! Read and understand these instructions before using the Blood Pressure Monitor. If you do not understand any part of these warnings, cautions or instructions, contact a healthcare professional for direction in the use of this product. If the Blood Pressure Monitor is not properly assembled and used, personal injury and damage to the Blood Pressure Monitor could result.**

⚠ **WARNING: DO NOT use this product without proper instruction from a healthcare professional.**

⚠ **WARNING: Operate the device only as intended. DO NOT use it for any other purpose.**

⚠ **WARNING: If components are damaged or missing, contact your Graham-Field authorized distributor immediately. DO NOT use substitute parts. Use only Lumiscope replacement parts. Non-Lumiscope replacement parts could cause personal injury and damage to the Blood Pressure Monitor.**


⚠ **WARNING: Contact your physician for specific information about your blood pressure. Self-diagnosis and treatment using measured results may be dangerous. Follow the instructions of your healthcare provider.**

⚠ **WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper assembly or use of this product.**

PRODUCT DESCRIPTION AND FEATURES

COMPONENTS INCLUDED WITH PRODUCT

1. Blood Pressure Monitor
2. Wrist Cuff
3. Arm Pillow / Storage Case
4. Gift Box
5. Instruction Manual

PRODUCT DESCRIPTION (SEE ABOVE)

1. Main Body
2. Display
3. Power Button (⏻)
4. Button M1
5. Button M2
6. Batteries (not included)
7. Arm Pillow / Storage Case
8. Battery cover
9. Wrist Cuff

CLASSIFICATION

- ME equipment not intended for use in an oxygen-rich environment or in the presence of flammable mixers
- Internally powered equipment
- Type BF equipment, applied part (cuff is recognized as applied part)

PRODUCT SYMBOLS

SYMBOL	MEANING
	Manufacturer
	Authorized Representative in the European community
	Symbol for the marking of electrical and electronics devices according to Directive 2002/96/EC. The device, accessories and packaging must be disposed of appropriately at the end of their use. Follow local ordinances and regulations for disposal.
	CE marking in conformity with EC directive 93/42/EEEC

	Keep dry
	Attention, consult accompanying documents
	Type BF Applied Part
	Power (ON/OFF) / Standby

SETUP – PREPARATION FOR OPERATION

BATTERY INSTALLATION / REPLACEMENT

- Replace the batteries when the battery indicator “” displays, or when nothing displays after the Power button is pressed.
- Use only fully-charged AAA alkaline batteries.
- DO NOT use rechargeable batteries.
- Always replace both batteries at the same time.
- Do not store the BPM with batteries if it will not be used for an extended period of time.

1. Remove the battery cover (Fig. 1).
2. Insert two AAA batteries into the battery compartment with polarities “+” and “-” matching polarity indicators.
3. Replace the battery cover.
4. Dispose of used batteries in compliance with local laws and regulations.

Fig. 1

SETTING THE DATE AND TIME

The Date and Time function provides an accurate time for each measurement. To get an accurate date and time, preset the date and time correctly before the first use of this device.

1. When the device is first operational, the display will show as in Fig. 2.
2. Press and hold the M1 button, then press the Power button, and the year number will flash as in Fig. 3.
3. Press M1 or M2 button to subtract or add the appropriate number, and press the Power button (⏻) to confirm.
4. When the year setup is finished, the month number will automatically flash. Press the M1 or M2 button to subtract or add the appropriate number, and press the Power button to confirm.
5. Press the Power button to finish setting the date and time. Follow the previous instructions to change the date and time.

Fig. 2

Fig. 3

Info: When in Date and Time mode, the device will automatically return to standby mode after one minute without operation.

PREPARATION FOR OPERATION

- Do not conduct any measurements if the temperature is low (below 50°F / 10°C) or high (over 104°F / 40°C), or if the relative humidity is higher than 85%, as this can lead to inaccurate readings.
- Take the measurement at room temperature in a quiet and stress-free environment.
- Do not move yourself or the unit during the measurement. Do not speak during the measurement.
- Blood pressure varies naturally depending on the time of day and is affected by many factors. Blood pressure is usually highest at work and reaches its lowest level during sleep.
- Blood pressure measurements should be assessed by a physician or trained healthcare professional who is familiar with your medical history. If you use the unit and regularly record the results, keep your physician informed of any changes in your blood pressure.
- The performance of this device can be affected as severe arrhythmias such as atrial or ventricular premature beats or atrial fibrillation are presented during measurement.
- The blood pressure measurements conducted with this unit are equivalent to measurements obtained by a trained observer in accordance with the values achieved using the cuff/stethoscope auscultation method and are within the specified EN 1060-4 standard limits.

OPERATION

OPERATION SAFETY

- ▲ **NOTICE:** The device is not waterproof. DO NOT immerse this device in liquid.
- ⚠ **WARNING:** Read and follow the entire instruction manual before operating this blood pressure monitor.
- ⚠ **WARNING:** DO NOT press the POWER button if the cuff has not been properly applied.
- ⚠ **WARNING:** DO NOT use the instrument if you think it is damaged or if you notice anything unusual.
- ⚠ **WARNING:** DO NOT disassemble or modify the device or cuff.
- ⚠ **WARNING:** DO NOT use this device if you have electrical implants.
- ⚠ **WARNING:** If you have had a mastectomy, do not use this device on the arm on the side of the mastectomy.
- ⚠ **WARNING:** DO NOT use this device on the wrist where intravascular access is present.
- ⚠ **WARNING:** DO NOT use this device simultaneously with other medical electrical equipment.
- ⚠ **WARNING:** DO NOT use this device in the presence of HF surgical equipment, MRI, or CT scanner.
- ⚠ **WARNING:** Discard old batteries carefully, out of reach of children. Swallowing the battery may be fatal. If the battery or other small parts are swallowed, contact a hospital immediately to have it removed.

- ⚠ **WARNING:** Avoid prolonged over-inflation of the bladder to prevent physical injury.
- ⚠ **WARNING:** If the cuff causes any discomfort during measurement, press the POWER button to turn off the device immediately.
- ⚠ **WARNING:** Keep blood pressure monitor away from children.
- ⚠ **WARNING:** This manual and the product are not substitutes for visiting the physician. Neither the information contained herein nor this product may be used to diagnose or treat health problems, or to prescribe drugs. If you have or suspect that you have a medical problem, seek immediate advice from your physician.
- ⚠ **WARNING:** Measuring too frequently may result in circulatory disorders, which can cause unpleasant sensations such as localized bleeding under the skin or temporary numbness in your arm. These symptoms do not usually last long; however, if you have not recovered quickly, consult your physician.
- ⚠ **WARNING:** Take into consideration the electromagnetic compatibility of the unit (e.g. disruptions to the power supply, radio frequency interference, etc.) Only use the unit indoors. To avoid inaccurate results due to electromagnetic interference between electrical and electronic equipment, do not use the unit near mobile phones or microwave ovens. Keep devices whose maximum power exceeds 2 W at least 11 feet (3.3 meters) from the blood pressure monitor.

- ▲ **NOTICE:** Avoid subjecting the monitor to shocks or vibrations, such as dropping it on the floor.

MEASUREMENT POSTURE

- Wait at least thirty minutes after finishing a caffeinated drink or cigarette before measurement.
- Sit down and relax for at least five minutes before measurement.
- Wait at least five minutes between repeated measurements.
- Use the same arm consistently for each measurement (preferably the left) and take the measurement at about the same time every day.
- Sit upright in a chair in a comfortable position with your feet flat on the floor in a natural position, your elbows placed on a table, and both feet on the ground. Do not interlock your legs during the measurement.
- Extend your measurement arm in front of you with the hand relaxed and the palm facing up.

APPLYING THE CUFF

- Ensure your wrist circumference is within applicable cuff range.


Fig. 4

Fig. 5

Fig. 6

- Remove all clothing from your arm, allowing the cuff to fit directly on the skin (fig. 4).

- Attach the cuff around your bare-skinned left wrist through the cuff with the display aligned with your palm. If it is not possible to take measurement with the left arm, use the right arm instead.
- Fasten the cuff with its bottom edge approximately 5/8 inch (~ 2 cm) above your palm line (Fig. 5). Ensure the cuff is not wrapped too tightly.
- Sit upright in a chair with your feet flat on the floor in a natural, comfortable position and relax. Rest your elbow on a table with the cuff at heart level and remain still (Fig. 6).

TAKING A MEASUREMENT

AUTOMATIC INFLATION

This device has four levels of inflation pressure: 190mmHg, 230mmHg, 270mmHg, and 300mmHg. When 190mmHg is insufficient, or wrist movement occurs, the device will automatically inflate to a higher pressure level to ensure successful measurement. **This is not a fault.**

RAPID DEFLATION DURING MEASUREMENT

If you do not feel well during measurement, or want to stop measurement for any reason, press the Power (⏻) button. The device will quickly release the air in the cuff and the device will return to standby mode.

MEASUREMENT


Fig. 7

Fig. 8

- Press the Power (⏻) button to start automatic measurement. The current date and time will display (Fig. 7).
- Press the Power button again; in two seconds all symbols will display, then "0 mmHg" will display. The pump will begin automatic inflation, and the inflation pressure, usually 190 mmHg, will display (Fig. 8).
- Press Button M1 or M2 to choose user number for memory storage. User selection can also be made by pressing Button M1 or M2 when readings are displayed after measurement is completed.

Fig. 9

Fig. 10


Fig. 11

- When the measurement is complete, the pump stops inflating and the pressure begins to decrease gradually, during which time the user's blood pressure and pulse will be calculated and display (Fig. 9).
- The cuff will then deflate quickly, and the blood pressure and pulse will display. M1 M2 will then flash to remind the user to record the reading (Fig. 10).
- Press button M1 or M2 to record the measurement in the appropriate memory. The measurement date and time will display alternately (Fig. 11). If M1 or M2 is not pressed, the measurement will not record.
- Press the POWER button to turn the device off, or it will turn off automatically after three minutes of non-operation. Readings will be saved into memory automatically while powering off.

MEASUREMENT INTERPRETATION

WHO Blood Pressure Classification

After each measurement is completed, the LCD will automatically display your position on the six segments of the bar indicator which correspond to the WHO (World Health Organization) Blood Pressure Indicator, shown above. Optimal and normal categories are indicated with 1 or 2 bars from the bottom which fall within the green indication; high-normal is indicated with the 3rd bar from the bottom which falls within the yellow indication; and grade 1-3 hypertension categories are indicated with 4-6 bars from the bottom which fall within the red indication. For example, if your blood pressure is 145mmHG / 88mmHg (systolic / diastolic), according to the WHO standard, your blood pressure level is Mild Hypertension.

Info: If the systolic and diastolic blood pressures fall into different categories, the higher value should be taken for classification. The WHO blood pressure classification indication is only a reminder; it can not be regarded as the final diagnosis. Always consult your physician for interpretation.

Irregular Heartbeat Detector

This device provides a blood pressure and pulse rate measurement even when an irregular heartbeat occurs. When the device detects an irregular heartbeat or any excessive body movement during measurement, the screen shown at upper right will display.

Info: Remain relaxed and still and do not talk during measurement

Info: Contact your physician if you see this indicator frequently.

Frequent and Low Battery Information

INDICATION	POSSIBLE REASON	CORRECTION METHOD
	The wrist cuff is incorrectly applied	Ensure the wrist cuff is correct applied and repeat measurement
	Arm or hand movement or talking occurred during measurement	Follow the recommendations in Taking a Measurement and repeat the measurement
	Weak batteries	Replace both batteries with new ones

IPD (Irregular Pulse Detection)

The device can detect irregular pulse (a pulse interval longer than 5/3 times the average pulse interval) during measurement. The IPD indicator will appear when more than three irregular pulses were detected during measurement. If the IPD indicator is displayed with measurement readings frequently, consult a physician for further direction.

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This function is equal to systolic minus diastolic. If systolic is higher than diastolic 60 mmHg, LCD will display "PP" and differentiation value. If PP indicator is displayed with measurement readings frequently, consult a physician for further direction.

MEMORY FUNCTION

Memory Recall


Fig. 12 Fig. 13 Fig. 14

- The device can store 90 sets of readings, respectively, in M1 and M2, and will automatically calculate the average value of the three latest readings for M1 and M2. When memory is full (90 sets of readings are stored), the oldest reading will be replaced by a new one. Readings in memory will not clear even if power is removed.
- When a measurement is finished, or the device is in standby mode, the user can press button M1 or M2 to recall memory. If you press button M1 or M2, the average value of the latest three readings (fig. 12) will display.
- Press again; "01" will display, the latest reading then change to another screen to display readings and measuring time as in fig. 13.
- Press again; "02" will display, which means second-to-last reading.

Memory Clearance

When a measurement is finished, or when the device is in standby mode, press and hold button M1 or M2 for at least five seconds; "CLR" will display, which means the stored reading for M1 or M2 is removed (fig. 14).

TROUBLESHOOTING

SYMPTOM	POSSIBLE CAUSE	REMEDY
No display after installing or replacing batteries	The batteries are depleted	Replace both batteries with new ones
	Battery polarity is incorrect	Install batteries with correct polarity
	Bad battery contact	Clean battery terminals with dry cloth

Inflation stops and re-inflates later	Automatic inflation should ensure correct measurement	See Automatic Inflation
	Arm or hand movement or talking during measurement	Remain still and silent during measurement
Reading is extremely low or high	Wrist cuff is not at same level as heart	Maintain correct posture during measurement
	Wrist cuff is improperly wrapped	Ensure wrist cuff is correctly wrapped
	Arm or hand strain during measurement	Remain still and silent during measurement
Pulse rate is too low or too high	Arm movement, hand movement, or talking occurred during measurement	Remain still and silent during measurement
	Batteries deplete quickly	Measurement taken immediately after exercise
	Faulty batteries	Use alkaline batteries of known manufacturers

MAINTENANCE

CLEANING

Clean the unit with a soft dry cloth.

▲ **NOTICE: DO NOT use any abrasive or volatile cleaners, solvents, naphtha, thinner, or gasoline to clean the device.**

▲ **NOTICE: DO NOT wet the cuff or attempt to clean it with water.**

▲ **NOTICE: DO NOT immerse any components of this device in liquid.**

STORAGE

Keep the device in its storage case when not in use.

▲ **NOTICE: Store the device in a clean, dry location, away from hot or cold temperature, humidity, or sunlight.**

▲ **NOTICE: Remove the batteries if the device will not be used for an extended period of time.**

SPECIFICATIONS

SPECIFICATION	DESCRIPTION		
Model	Lumiscope 1143		
Dimensions (LxWxH)	2.17 in. x 2.76 in. x 1.18 in. / 66 mm x 70 mm x 30 mm		
Weight (without batteries)	~ 3.88 oz. / 110 g		
Measuring Method	Oscillometric		
Extreme Pressure / cuff Pressure	300mmHg		
Measuring Range	DIA	40mmHg - 180 mmHg	
	SYS	60 - 260 mmHg	
	Pulse Rate	40 - 160 beats per minute	
Measuring Accuracy	Static Pressure	±3mmHG	
	Pulse Rate	±5% reading	
Inflation	Automatic by pump		
Rapid Deflation	Automatic electronic valve		
Batteries	Two AAA x 1.5V (not included)		
Environment	Operation	Temperature	10°C - 40°C
		Humidity	≤ 85%
		Air Pressure	700hPa - 1060hPa
	Transport and Storage	Temperature	-20°C - 50°C
		Humidity	≤ 85%
		Air Pressure	500hPa - 1060hPa
Cuff Size	Appropriate for wrist circumference 4.92 in. ~ 8.07 in. / 12.5 cm ~ 20.5 cm		
Complete Kit	BPM body, Wrist cuff, Instruction manual, Gift box, and Storage case		
Pollution Degrees	2 degrees		
Overvoltage category	Category II		
High Altitude	≤ 6,562 ft. / 2,000 m		

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a component is warranted by a third party, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported to GF's customer service team within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. This limited warranty is not transferable. Within the guidelines set forth in this document, the following components are warranted for the time period set forth below:

COMPONENTS

Digital Blood Pressure Monitor.....24 months
Cuff.....1 year

* Labor is not included in the warranty

† Upholstery is only warranted on material supplied by GF: 1 year

‡ The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified.

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

This limited warranty shall only apply to defects that are reported to the Distributor from whom the Customer purchased the product within the applicable warranty period. If there is not a Distributor, you must contact Hausted directly by calling 1.678.291.3207, or by e-mailing a request to cs@grahamfield.com. Specific directions will be provided by the Customer Service Representative. Failure to abide by the specific directions will result in denial of the warranty claim.

EXCLUSIONS

The warranty does not cover and Hausted shall not be liable for the following:

- Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- Products considered to be of a non-durable nature including, but not limited to: tips, casters, filters, fuses, gaskets, lubricants, and charts;
- Accessories or parts not provided by Intensa;
- Matching of color, grain or texture except to commercially acceptable standards;
- Changes in color caused by natural or artificial light;
- Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by Intensa;
- Any labor or shipping charges incurred in the replacement part installation or repair;
- Costs and expenses of regular maintenance and cleaning; and
- Representations and warranties made by any person or entity other than Intensa.

The warranty does not cover and GF shall not be liable for the following:

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS.

The warranties contained herein contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document. For additional information on this product or this warranty, contact a GF Customer Service Representative.

NOTES:

- Additional terms and conditions may apply.
- Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- Claims for any short shipment must be made within three (3) days of the invoice date.

LUMISCOPE
by graham-field
grahamfield
Manufactured For:
GF Health Products, Inc.
One Graham-Field Way, Atlanta, GA 30340

*1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website:
www.grahamfield.com

© 2009, GF Health Products, Inc. All Rights Reserved.
Lumiscope is a registered trademark of GF Health Products, Inc.

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© 2009, GF Health Products, Inc. Todos los derechos reservados. Lumiscope es una marca de GF Health Products, Inc.