

- ⚠ **WARNING: Important! Read and understand these instructions before assembling or using the UpRise Onyx Folding Walker. If you do not understand any part of these warnings or instructions, contact a healthcare professional for direction in the use of this product. If the UpRise Onyx Folding Walker is not properly assembled, personal injury and damage to the Walker could result.**
- ⚠ **WARNING: If components are damaged or missing, contact your dealer immediately. DO NOT use substitute parts. Use only Lumex replacement parts. The use of non-Lumex replacement parts could cause personal injury, property damage, and void the warranty.**
- ⚠ **WARNING: This product has a maximum weight capacity of 400 lb (180 kg), EVENLY DISTRIBUTED.**
- ⚠ **WARNING: Always use caution when using this product on wet or slippery surfaces.**
- ⚠ **WARNING: GF Health Products, Inc. ("GF") assumes no responsibility for any damage or injury caused by improper assembly or use of this product.**

INTENDED USE

The Lumex 700175CR UpRise Onyx Folding Walker is intended as a dual-purpose mobility aid and rising aid. It is recommended for use by users whose height is between 5'4" and 6'2". Please consult a healthcare professional for proper adjustment, as recommended user heights are an approximation.

UNPACKING

Remove walker frame and leg attachments from carton and protective packaging. Check for obvious damage to the carton or its contents. If damage is evident, do not use the product; notify the carrier and your Graham-Field authorized distributor for further instruction.

ASSEMBLY

Install Leg Attachments

1. Depress the leg assembly button and insert the leg attachment into the leg until the leg assembly button pops out of its hole on the frame, as shown at right.
2. Repeat for other three legs.

Adjust Walker Height

Handgrip-to-floor height is adjustable between 32" and 39". Each walker leg is equipped with several height adjustment positions, as shown at right.

1. To adjust leg height, depress the height adjustment button and twist the leg attachment to the desired height until the height adjustment button pops out of the appropriate hole.
2. Repeat for other three legs.


- ⚠ **WARNING: Before proceeding, ensure leg assembly buttons protrude fully from frame holes, height adjustment buttons protrude fully from height adjustment holes, and all legs are adjusted to the same height.**

Open Walker (see picture at right)

1. Grasp one handgrip and pull away from center of walker.
2. Repeat with other handgrip.
3. Push handgrips away from each other **until you hear the walker release mechanism click into place.**
4. Walker is properly opened when release mechanism push button is up.

Close Walker

1. Push downward on release mechanism and hold—sides of walker will collapse inward when release mechanism is depressed.
2. Grasp one handgrip and push toward center of walker.
3. Grasp the other handgrip and push toward center of walker.



700175CR Lumex UpRise Onyx Folding Walker

WALKER OPERATION

⚠ WARNING: Before use, ensure walker is fully opened, leg assembly buttons protrude fully from frame holes, height adjustment buttons protrude fully from height adjustment holes, and all legs are adjusted to the same height.

Once leg attachments are adjusted and the walker is opened, the user can operate the walker by grasping one hand grip in each hand, positioning body between the two rear walker legs, and ambulating. The rear walker legs should be in contact with the floor at all times.

PORTABLE TOILET SAFETY FRAME OPERATION

The walker may also be used as a portable toilet safety frame in either position shown at right.



MAINTENANCE

⚠ WARNING: Inspect the walker weekly for proper function and wear. If any component is not in proper working order or appears worn, contact your dealer immediately for repair or appropriate replacement parts.

Wipe walker weekly with a clean, damp cloth. Clean handgrips with mild soap or detergent.

▲ NOTICE: DO NOT use solvents, abrasive cleaners, or scouring pads on any part of the walker.

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. (See Obtaining Warranty Service below) This limited warranty is not transferable.

The warranted components and time periods are set forth below:

Frame:three (3) years
Normal wear components such as tips, wheels/casters, and handgrips:three (3) months

Labor is not included in the warranty.

The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified.

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1.770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim.

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Matching of color, grain or texture except to commercially acceptable standards;
- 6) Changes in color caused by natural or artificial light;
- 7) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 8) Any labor or shipping charges incurred in the replacement part installation or repair;
- 9) Costs and expenses of regular maintenance and cleaning; and
- 10) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.

LUMEX[®]
by graham-field

Manufactured For:

GF Health Products, Inc.
One Graham-Field Way, Atlanta, GA 30340
Made in China

© 2009, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

LX_700175CR-INS-LAB-RevC23



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

graham
field