Transcutaneous Electrical Nerve Stimulation (TENS)

Model GF-3 / GF-3T

Operation Manual

Read this manual before operating your GF-3 / GF-3T.

Save this manual for future use.

CONTENTS

GENERAL DESCRIPTION	3
WHAT IS TENS?.....	3
INDICATIONS AND CONTRAINDICATIONS.....	3
SAFETY	4
ABOUT THE DEVICE.....	6
EXPLANATION OF CONTROL FUNCTIONS.....	7
ATTACHING THE LEAD WIRES	8
ELECTRODE SELECTION AND CARE	8
TIPS FOR SKIN CARE	9
CONNECTING THE DEVICE.....	10
BATTERY INFORMATION	11
CARING FOR YOUR DEVICE.....	12
TROUBLESHOOTING.....	13
TECHNICAL SPECIFICATIONS	14
OUTPUT SPECIFICATIONS	14
LIMITED WARRANTY	15
MANUEL D'UTILISATION EN FRANÇAIS	17
MANUAL DE OPERACIÓN EN ESPAÑOL.....	33

GENERAL DESCRIPTION

TENS, Transcutaneous Electrical Nerve Stimulation, is a method of relieving symptomatic chronic intractable pain.

WHAT IS TENS?

TENS is a treatment whereby electrical impulses are applied to nerves through electrode pads placed on the skin. TENS is non-invasive and does not use pharmaceuticals.

TENS uses a two-pronged approach to pain relief. First, sensory nerves are targeted, stimulating them to block pain signals and prevent their transmission to the brain. Second, TENS promotes the production of endorphins — neurochemicals occurring naturally in the brain — which have analgesic properties.

INDICATIONS AND CONTRAINDICATIONS

Read the operation manual before using this TENS device.

Federal law (USA) restricts this device to sale by or on the order of a physician.

Indications


TENS is indicated to be used under a physician's prescription for the symptomatic relief of chronic intractable pain.

Contraindications




- Any electrode placement which applies current to the carotid (neck) region.
- Patients with implanted electronic devices (for example, a pacemaker) or metallic implants should not undergo TENS treatment without first consulting a physician.
- Any electrode placement which causes current to flow transcerebrally (through the head).
- The use of unit whenever pain symptoms are undiagnosed and the etiology is unknown.

SAFETY

Always follow basic safety precautions, including the following:

-  **WARNING:** Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious personal injury.
-  **CAUTION:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in moderate or minor personal injury.
-  **NOTICE:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in product or property damage.

WARNINGS

-  **WARNING: Explosion hazard**
Explosion hazard is possible if used in the presence of explosives, flammable materials or flammable anesthetics.
-  **WARNING: Heart disease**
Caution should be used when applying the device to patients suspected of having heart disease. Further clinical data is needed to show if there are adverse side effects on individuals with heart disease.
-  **WARNING: Keep this device out of the reach of children.**
-  **WARNING: The safety of the device during pregnancy or delivery has not been established.**
-  **WARNING: Do not place electrodes on front of the throat. This may result in spasms of the laryngeal and pharyngeal muscles.**
-  **WARNING: Do not place the electrodes over the carotid nerve.**
-  **WARNING: The device is not effective for pain of central origin (headaches).**
-  **WARNING: Avoid adjusting controls while operating machinery or vehicles.**
-  **WARNING: The device may interfere with electronic monitoring equipment (such as ECG monitors and ECG alarms).**

- ⚠ WARNING: Do not change any mode during treatment.**
- ⚠ WARNING: Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. Using an alternate conductive medium or alternate electrode placement can usually reduce the irritation. Consult your physician / clinician before using an alternative conductive medium or electrode placement.**
- ⚠ WARNING: Electrodes should not be placed over the eyes, in the mouth, or internally.**
- ⚠ WARNING: The device has no curative value.**
- ⚠ WARNING: TENS devices should be used only under the continued supervision of a physician / clinician.**
- ⚠ WARNING: TENS is a symptomatic treatment and as such suppresses the sensation of pain which would otherwise serve as a protective mechanism.**

Cautions and Adverse Reactions

- ⚠ CAUTION: Isolated cases of skin irritation may occur at the site of electrode placement following long-term application.**
- ⚠ CAUTION: If skin irritation occurs TENS treatment should be stopped and electrodes removed until the cause of the irritation can be determined.**
- ⚠ CAUTION: Effectiveness is highly dependent upon patient selection of a doctor qualified in the management of pain patients.**
- ⚠ CAUTION: If the device treatment becomes ineffective or unpleasant, stimulation should be discontinued until reevaluation by a physician / clinician.**
- ⚠ CAUTION: Always turn the device OFF before applying or removing electrodes.**

ABOUT THE DEVICE

This device is a battery-operated device that includes two controllable output channels. This device creates electrical impulses in which intensity, duration, and modulation can be altered. The device controls are easy to use and the slide cover protects accidental changes in settings.


System Components

Your device will include the following components or accessories:

- TENS unit
- Carrying case
- Lead wires
- Operation manual
- Electrodes

Also required (not included): One 9-Volt alkaline battery.

Device Controls

Front Panel Cover

This cover located on the front of the unit conceals the controls for Pulse Width, Pulse Rate, Mode Selector, and Timer (GF-3T only). Press the top portion of the cover and pull down in order to open the cover.

EXPLANATION OF CONTROL FUNCTIONS

Ch1 / Ch2 Knobs	The Channel 1 and Channel 2 Intensity Control Knobs control the strength of the stimulation and also function as ON/OFF controls for each channel.	
MODE Selection Switch	The MODE Selection Switch sets the treatment mode. The three modes are <u>Burst</u> (B), <u>Continuous</u> (C), and <u>Modulation</u> (M).	<u>Burst Mode</u> (B) releases individual bursts twice per second, pulse width is adjustable, and the pulse rate is set at 100 Hz per second.
		<u>Continuous Mode</u> (C) delivers stimulation continuously at the settings determined by the Intensity, Width, and Rate Knobs.
		<u>Modulation Mode</u> (M) decreases the pulse width to 60% of the original setting; this decreased pulse width is maintained for 1.5 seconds before returning to the original pulse width setting, which is maintained for 3.5 seconds. The cycle is then repeated. The intensity and pulse rate are adjustable.
Pulse Width Control Knob	The Pulse Width Control Knob regulates the pulse width for both channels.	
Pulse Rate Control Knob	The Pulse Rate Control Knob regulates the number of pulses per second for both channels.	
Timer Selection Switch * (GF-3T Only)	The Timer Selection Switch selects the duration of TENS treatment. The three Timer selections are <u>Continuous</u> (C), <u>30 minutes</u> (30), and <u>60 minutes</u> (60). * To resume operation or reset the Timer, turn the Intensity Control Knob OFF and then ON.	

ATTACHING THE LEAD WIRES

- ⚠ WARNING:** Ensure the device is OFF before connecting the lead wires.
- ⚠ WARNING:** Never insert the lead wire plug into an AC power supply socket. Personal injury and/or damage to the TENS unit could occur.
- ▲ NOTICE:** Use care when you plug and unplug the wires. Pulling on the lead wire instead of its insulated connector may cause wire breakage.

The lead wires provided with the device insert into the ports located the top of the unit. After connecting the wires to the unit, attach each wire to an electrode.

Lead wires provided with the device are compliant with mandatory compliance standards set forth by the FDA.

ELECTRODE SELECTION AND CARE

Using Electrodes

Use the electrodes as prescribed. Follow application procedures outlined in electrode packaging to maintain stimulation and prevent skin irritation. The electrode packaging provides instructions for care, maintenance, and proper storage of electrodes.

TIPS FOR SKIN CARE

Good skin preparation is important for effective and comfortable use of your TENS device.

- Always clean the electrode site with mild soap and water solution, rinse well, and dry thoroughly prior to any electrode application.
- Any excess hair should be clipped, not shaved, to ensure good electrode contact with the skin.
- If a skin treatment or preparation is recommended by your physician / clinician, apply the skin treatment as recommended, let dry, and apply electrodes as directed. Following these recommendations will both reduce the chance of skin irritation and extend the life of your electrodes.
- Avoid excessive stretching of the skin when applying electrodes. Proper application is best accomplished by applying the electrode, then smoothly pressing it in place from the center outward.
- When removing electrodes, always remove by pulling in the direction of hair growth.
- It may be helpful to rub skin lotion on electrode placement area when not wearing electrodes.

CONNECTING THE DEVICE

Insert battery

Turn the device to the OFF position before inserting or removing the battery. When inserting the battery, ensure the battery polarity (+ and -) markings match the markings on the device.

Prepare the Skin

Prepare the skin as previously described and according to the instructions provided with your electrodes. Before attaching the electrodes, identify the area that your physician / clinician has recommended for electrode placement.

1. Connect the lead wires to the electrodes: connect the lead wires to the electrodes before applying the electrodes to the skin.

WARNING: Ensure both Intensity Control Knobs for Channel 1 and 2 are turned to the OFF Position (counterclockwise) before applying the electrodes.

2. Place electrodes on the skin: place the electrodes on the skin as recommended by your physician / clinician.
3. Insert lead wire connector into the device: plug end of lead wire into the channel output port (jack) to be used; push the plug in as far as it will go.
4. Select treatment settings: ensure your unit is still set to the proper settings recommended by your physician / clinician.
5. Adjusting Channel Intensity Control: locate the Intensity Control Knob (Channel 1 or 2) at the top of the unit. Slowly turn the Intensity Control Knob clockwise until the stimulation is at the level recommended by your physician / clinician (if you don't feel anything, turn the Knob OFF then ON again and carefully turn the Control Knob until you feel a tingling or slight twitch under or around the electrodes). Always start with the lowest setting and increase the intensity slowly.

If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level; or cease stimulation and contact your physician.

BATTERY INFORMATION

When the indicator lights on the front of the unit no longer illuminate, the battery has become too weak to power the unit, and the existing battery should be replaced with a new battery. At this point, the unit will turn OFF until a new battery is inserted.

▲ NOTICE: GF Health Products, Inc. recommends the use of only a 9V alkaline battery with this device.

Replacing the Battery

When the indicator lights on the front of the unit no longer illuminate, the battery should be replaced.

1. Turn unit OFF.
2. Remove the front panel cover by pressing on the top of the panel and pressing down in order to slide the panel down. Continue sliding the panel downwards until the panel is completely removed from the unit. This will reveal the battery compartment.
3. Remove the discharged battery from the device.
4. Place new battery in the compartment. Note: Be sure the proper polarity (+ and -) markings match the markings in the device.
5. Replace the front panel cover.
6. Dispose of the old battery according to local guidelines and regulations.

CARING FOR YOUR DEVICE

Your device may be cleaned by wiping gently with a damp cloth moistened with mild soap and water. Do not immerse the device in water or other liquids.

Wipe lead wires with damp cloth moistened with soap and water. Do not immerse the lead wires.

To properly store the device for an extended period of time, remove the battery from the unit. Place the unit and accessories in the carrying case provided and store in a cool, dry location.

TROUBLESHOOTING

If the device does not function properly:

1. Ensure the battery is properly installed or replace the battery. Be sure to observe proper polarity markings when replacing the battery. If the indicator lights on the front of the unit no longer illuminate when it is turned on, replace the battery and check again.
2. If the indicator light is flashing, the intensity has been adjusted, and no stimulation is felt, check to ensure the lead wires are properly connected and the electrodes are properly applied to the skin. If the unit appears to be functioning and no stimulation is felt, the lead wires or electrodes may need to be replaced.
3. If the battery appears to be charged and the unit is not functioning, turn both Intensity Control Knobs to the OFF position (counterclockwise). Then gradually turn the Intensity Control Knob (clockwise) until stimulation is felt. If device still is not working, turn the unit off and contact your authorized GF Health Products, Inc. distributor.

If there is any other problem, please contact an authorized GF Health Products, Inc. distributor. Do not try to repair a defective device.

TECHNICAL SPECIFICATIONS

Channel	Dual, isolated between channels
Modes of Operation	Burst, Continuous, and Modulation
Pulse Intensity	Adjustable 0 mA – 80 mA peak into 500 ohm load each channel, constant current
Pulse Width	30uS-260uS (adjustable)
Pulse Rate	2 Hz-150 Hz (adjustable)
Timer (GF-3T Only)	Continuous, 30 minutes, and 60 minutes
Burst Mode (B)	Burst consists of 2 bursts per second at 100 Hz
Wave Form	Asymmetrical bi-phasic square pulse
Voltage	0 – 100 Volt (open circuit)
Power Source	9-Volt battery
Dimensions	95 mm x 65 mm x 23.5 (height x width x thickness)
Weight	115 g (including battery)

OUTPUT SPECIFICATIONS

Mode	Intensity (mA)	Pulse Width (uSec)	Pulse Rate Frequency (Hz)	Cycle Time (Sec)
Burst	Adj. 0-80	Adj. 30-260	100 Hz fixed N/A 2 bursts per Sec	
Continuous	Adj. 0-80	Adj. 30-260	Adj. 2-150 Hz	N/A
Modulation	Adj. 0-80	Modulates down from preset width setting by 60% then back to original setting	Adj. 2-150 Hz	5 Sec total time

**LIMITED WARRANTY:
SCOPE OF WARRANTY**

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. (See Obtaining Warranty Service below) This limited warranty is not transferable.

The warranted components and time periods are set forth below:

TENS Unit: 1 (one) year

* Labor is not included in the warranty.

‡ The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified. The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1.770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim.

The warranty does not cover and GF shall not be liable for the following:

1. Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
3. Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
4. Accessories or parts not provided by GF;
5. Matching of color, grain or texture except to commercially acceptable standards;
6. Changes in color caused by natural or artificial light;
7. Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
8. Any labor or shipping charges incurred in the replacement part installation or repair;
9. Costs and expenses of regular maintenance and cleaning; and
10. Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document. For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

1. Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com.
2. Freight claims must be noted on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
3. Claims for any short shipment must be made within three (3) days of the invoice date.

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© 2008, GF Health Products, Inc. Tous droits réservés.