

Installation tools are provided (open-end M10 x M11 wrench and M4 hex key)

- ⚠ **WARNING: Important! Read and understand these instructions before assembling or using the 8146 Utility Cart. If you do not understand any part of these instructions, contact your GF Health Products, Inc. ("GF") authorized distributor for direction in the use of this product. If the Utility Cart is not properly assembled, personal injury and/or damage to the Utility Cart could result.**
- ⚠ **WARNING: If components are damaged or missing, contact your GF authorized distributor immediately. DO NOT use substitute parts. Use only Graham-Field replacement parts.**
- ⚠ **WARNING: The 8146 Utility Cart has a maximum load capacity of 400 lb (181 kg).**
- ⚠ **WARNING: GF Health Products, Inc. assumes no responsibility for any damage or injury caused by improper assembly or use of this product.**

ASSEMBLY (REFER TO PICTURES AT RIGHT)

- ⚠ **WARNING: Assembly, repair and maintenance must be performed by qualified personnel only.**

The 8146 Utility Cart ships unassembled, as shown at right. Included components are listed below. Items 3, 4, 5, 6 and 7 are packaged in a bagged hardware kit.

8146 Utility Cart Components

Item	Qty	Description	Item	Qty	Description
1	3	Shelf	4	12	Lock Nut M6
2	2	Leg (with Casters)	5	12	Lock Washer M6
3	12	Button Socket Head Cap Screw M6 x 35	6	1	Open-End Wrench
			7	1	M4 Hex Key

1. Carefully remove packing from components.
2. Place shelves between legs, ensuring that handles at tops of legs face outward, and that shelves are positioned as shown at lower right.
3. Align leg holes with shelf holes as shown.
4. Install lock washers on screws. From outside, insert screws through leg holes and shelf holes. Install lock nuts on screws.
5. Ensure the legs are attached to the shelves in an upright perpendicular position and the shelves are level as shown.
6. Use the provided tools to tighten the screws securely.

OPERATION


- ▲ **NOTICE: Ensure all components are secure before each use.**

The Utility Cart may be pushed from either end.

MAINTENANCE, CLEANING AND STORAGE

- ▲ **NOTICE: Replace cracked, damaged, corroded or damaged components immediately.**
- ▲ **NOTICE: Do not use acid, aggressive disinfectant or corrosive cleaner to clean this product.**
- ▲ **NOTICE: Ensure Utility Cart is completely dry before packing or storing, and store only in dry conditions.**

To clean Utility Cart, wipe with a clean, soft paper or cloth. To disinfect Utility Cart, wipe with a clean, soft, paper or cloth dampened with a general purpose, noncorrosive cleaner, rinse with clear water, and dry with a clean, soft paper or cloth.



8146 Utility Cart unassembled



8146 Utility Cart assembled

WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted. This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. (See Obtaining Warranty Service below) This limited warranty is not transferable.

The warranted components and time periods are set forth below:

Utility Cart..... 1 (one) year

Labor is not included in the warranty.

The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified.

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1.770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim.

The warranty does not cover and GF shall not be liable for the following:

- 1) Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- 2) Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
- 3) Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
- 4) Accessories or parts not provided by GF;
- 5) Matching of color, grain or texture except to commercially acceptable standards;
- 6) Changes in color caused by natural or artificial light;
- 7) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
- 8) Any labor or shipping charges incurred in the replacement part installation or repair;
- 9) Costs and expenses of regular maintenance and cleaning; and
- 10) Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document.

For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

- 1) Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com.
- 2) Freight claims must be notated on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
- 3) Claims for any short shipment must be made within three (3) days of the invoice date.



Manufactured For:

GF Health Products, Inc.
1 Graham Field Way, Atlanta, GA 30340
Made in India



+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.
www.grahamfield.com

© 2008, GF Health Products, Inc. All Rights Reserved.
Graham-Field is a trademark of GF Health Products, Inc.

GO_8146-INS-LAB-RevE25

