Transcutaneous Electrical Nerve Stimulation (TENS), Digital Model GF-TX5

Operation Manual

Read this manual before operating your GF-TX5.

Save this manual for future use.

1639

CONTENTS

GENERAL DESCRIPTION.....	3
WHAT IS TENS?	3
INDICATIONS AND CONTRAINDICATIONS.....	3
SAFETY	4
ABOUT THE DEVICE	6
EXPLANATION OF KEY / KNOB CONTROL FUNCTIONS	7
ATTACHING THE LEAD WIRES	8
ELECTRODE SELECTION AND CARE	8
TIPS FOR SKIN CARE	8
CONNECTING THE DEVICE.....	9
BATTERY INFORMATION.....	10
CARING FOR YOUR DEVICE	11
TROUBLESHOOTING	11
TECHNICAL SPECIFICATIONS.....	12
LIMITED WARRANTY	15
MANUEL D'UTILISATION EN FRANÇAIS.....	17
MANUAL DE OPERACIÓN EN ESPAÑOL	33

GENERAL DESCRIPTION

TENS, Transcutaneous Electrical Nerve Stimulation, is a method of relieving symptomatic chronic intractable pain.

This unit is a dual-channel digital stimulator for active treatment application, which has a Liquid Crystal Display indicating operation modes and output as well as an 8-bit microcomputer for controlling the system.

The electronics of the unit create electric impulses; the intensity, duration, frequency per second and modulation of these impulses can be adjusted.

WHAT IS TENS?

TENS is a treatment whereby electrical impulses are applied to nerves through electrode pads placed on the skin. TENS is non-invasive and does not use pharmaceuticals.

TENS uses a two-pronged approach to pain relief. First, sensory nerves are targeted, stimulating them to block pain signals and prevent their transmission to the brain. Second, TENS promotes the production of endorphins — neurochemicals occurring naturally in the brain — which have analgesic properties.

INDICATIONS AND CONTRAINDICATIONS

Read the operation manual before using this TENS device.

Federal law (USA) restricts this device to sale by or on the order of a physician.

Indications

TENS is indicated to be used under a physician's prescription for the symptomatic relief of chronic intractable pain.

Contraindications

- Any electrode placement which applies current to the carotid (neck) region.
- Patients with implanted electronic devices (for example, a pacemaker) or metallic implants should not undergo TENS treatment without first consulting a physician.
- Any electrode placement which causes current to flow transcerebrally (through the head).
- The use of unit whenever pain symptoms are undiagnosed and the etiology is unknown.

SAFETY

Always follow basic safety precautions, including the following:

-  **WARNING:** Indicates a potential hazard situation or unsafe practice that, if not avoided, could result in death or serious personal injury.
-  **CAUTION:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in moderate or minor personal injury.
-  **NOTICE:** Indicates a potential hazard or unsafe practice that, if not avoided, could result in product or property damage.

Warnings

-  **WARNING: Explosion hazard:** Explosion hazard is possible if used in the presence of explosives, flammable materials or flammable anesthetics.
-  **WARNING: Heart disease:** Caution should be used when applying the device to patients suspected of having heart disease. Further clinical data is needed to show if there are adverse side effects on individuals with heart disease.
-  **WARNING:** Keep this device out of the reach of children.
-  **WARNING:** The safety of the device during pregnancy or delivery has not been established.
-  **WARNING:** Do not place electrodes on front of the throat. This may result in spasms of the laryngeal and pharyngeal muscles.
-  **WARNING:** Do not place the electrodes over the carotid nerve.
-  **WARNING:** The device is not effective for pain of central origin (headaches).
-  **WARNING:** Avoid adjusting controls while operating machinery or vehicles.
-  **WARNING:** The device may interfere with electronic monitoring equipment (such as ECG monitors and ECG alarms).

-  **WARNING:** Do not change any mode during treatment.
-  **WARNING:** Some patients may experience skin irritation or hypersensitivity due to the electrical stimulation or electrical conductive medium. Using an alternate conductive medium or alternate electrode placement can usually reduce the irritation. Consult your physician / clinician before using an alternative conductive medium or electrode placement.
-  **WARNING:** Electrodes should not be placed over the eyes, in the mouth, or internally.
-  **WARNING:** The device has no curative value.
-  **WARNING:** TENS devices should be used only under the continued supervision of a physician / clinician.
-  **WARNING:** TENS is a symptomatic treatment and as such suppresses the sensation of pain which would otherwise serve as a protective mechanism.

Cautions and Adverse Reactions

-  **CAUTION:** Isolated cases of skin irritation may occur at the site of electrode placement following long-term application.
-  **CAUTION:** If skin irritation occurs TENS treatment should be stopped and electrodes removed until the cause of the irritation can be determined.
-  **CAUTION:** Effectiveness is highly dependent upon patient selection of a doctor qualified in the management of pain patients.
-  **CAUTION:** If the device treatment becomes ineffective or unpleasant, stimulation should be discontinued until reevaluation by a physician / clinician.
-  **CAUTION:** Always turn the device OFF before applying or removing electrodes.

ABOUT THE DEVICE

This device is a battery-operated device that includes two controllable output channels. This device creates electrical impulses in which intensity, duration, and modulation can be altered. The device controls are easy to use and the slide cover protects accidental changes in settings.


System Components

Your device will include the following components or accessories:

- TENS unit
- Carrying case
- Lead wires
- Operation manual
- Electrodes

Also required (not included): One 9-Volt alkaline or nickel-cadmium rechargeable battery.

Device Controls

Slide Cover

This cover located on the front of the unit conceals the controls for DOWN, MODE, SET, and UP. Press the top portion of the cover and pull down in order to open the cover.

EXPLANATION OF KEY / KNOB CONTROL FUNCTIONS

DOWN Key	This key decreases the setting when in <u>Pulse Width Mode</u> or <u>Pulse Rate Mode</u> . *	This key regulates the number of pulse width or pulse rate of the individual current pulses.
	*Decrease by pressing the DOWN key: The width can be adjusted by 10uS/step. The rate can be adjusted in 1-20Hz by 1Hz/step, 20Hz-150Hz by 5Hz/step.	
MODE Key	This key selects <u>Timer Mode</u> , <u>Stimulation Mode</u> , <u>Alternate Mode</u> , <u>Pulse Width Mode</u> (uS) or <u>Pulse Rate Mode</u> (Hz).	This key changes the treatment parameter. Each time the MODE key is pressed, the next parameter will display. The selected treatment parameter in the current mode will flash.
SET Key	This key switches between the different settings in the <u>Timer Mode</u> , <u>Stimulation Mode</u> , and <u>Alternate Mode</u> treatment parameters.	Each time the SET key is pressed, the parameter will change to the next setting. The selected mode setting will flash. When the desired parameter is flashing, press the MODE key to switch to the next treatment parameter. The parameter just set will be displayed and will no longer flash.
UP Key	This key increases the setting when in <u>Pulse Width Mode</u> or <u>Pulse Rate Mode</u> . *	This key regulates the number of pulse width or pulse rate of the individual current pulses.
	*Increase by pressing the UP key: The width can be adjusted by 10uS/step. The rate can be adjusted in 1-20Hz by 1Hz/step, 20Hz-150Hz by 5Hz/step.	
Ch1 / Ch2 Knobs	Channel 1 and Channel 2 Intensity Control Knobs	These knobs control the strength of the stimulation and also function as ON / OFF controls.

ATTACHING THE LEAD WIRES

 **WARNING:** Ensure the device is OFF before connecting the lead wires.

 **WARNING:** Never insert the lead wire plug into an AC power supply socket. Personal injury and/or damage to the TENS unit could occur.

▲ NOTICE: Use care when you plug and unplug the wires. Pulling on the lead wire instead of its insulated connector may cause wire breakage.

The lead wires provided with the device insert into the ports located on the top of the unit. If only one lead will be used, plug it into the channel 1 port. After connecting the wires to the unit, attach each wire to an electrode.

Lead wires provided with the device are compliant with mandatory compliance standards set forth by the FDA.

ELECTRODE SELECTION AND CARE

Using Electrodes

Use the electrodes as prescribed. Follow application procedures outlined in electrode packaging to maintain stimulation and prevent skin irritation. The electrode packaging provides instructions for care, maintenance, and proper storage of electrodes.

TIPS FOR SKIN CARE

Good skin preparation is important for effective and comfortable use of your TENS device.

- Always clean the electrode site with mild soap and water solution, rinse well, and dry thoroughly prior to any electrode application.
- Any excess hair should be clipped, not shaved, to ensure good electrode contact with the skin.
- If a skin treatment or preparation is recommended by your physician / clinician, apply the skin treatment as recommended, let dry, and apply electrodes as directed. Following these recommendations will both reduce the chance of skin irritation and extend the life of your electrodes.
- Avoid excessive stretching of the skin when applying electrodes. Proper application is best accomplished by applying the electrode, then smoothly pressing it in place from the center outward.

- When removing electrodes, always remove by pulling in the direction of hair growth.
- It may be helpful to rub skin lotion on electrode placement area when not wearing electrodes.

CONNECTING THE DEVICE

Inserting the Battery

Turn the device to the OFF position before inserting or removing the battery. When inserting the battery, ensure the battery polarity (+ and -) markings match the markings on the device.

Preparing the Skin

Prepare the skin as previously described and according to the instructions provided with your electrodes. Before attaching the electrodes, identify the area that your physician / clinician has recommended for electrode placement.

1. **Connect the lead wires to the electrodes:** connect the lead wires to the electrodes before applying the electrodes to the skin.

WARNING: Ensure both Intensity Controls for Channel 1 and 2 are turned to the OFF position (counterclockwise) before applying the electrodes.

2. **Place electrodes on the skin:** place the electrodes on the skin as recommended by your physician / clinician.
3. **Insert lead wire connector into the device:** plug end of lead wire into the channel output port (jack) to be used; push the plug in as far as it will go.
4. **Select treatment settings:** ensure your unit is still set to the proper settings recommended by your physician / clinician.
5. **Adjusting Channel Intensity Control:** locate the Intensity Control Knob (Channel 1 or 2) at the top of the unit. Slowly turn the Intensity Control Knob clockwise until the stimulation is at the level recommended by your physician / clinician. (If you don't feel anything, turn the Knob OFF then ON again and carefully turn the Control Knob until you feel a tingling or slight twitch under or around the electrodes.) Always start with the lowest setting and increase the intensity slowly.

If the stimulation levels are uncomfortable or become uncomfortable, reduce the stimulation intensity to a comfortable level; or cease stimulation and contact your physician.

6. **Setting the Patient Compliance Counter:**

- a) To turn on the Patient Compliance Counter: While the unit is ON, hold down the UP button and press the MODE button simultaneously.
- b) To reset the Patient Compliance Counter: While the unit is ON, press the UP button and the MODE button simultaneously (this will take you into the Compliance Counter), then push the DOWN button and press the MODE button simultaneously.
- c) Press the UP and MODE buttons simultaneously to return to the treatment status.

BATTERY INFORMATION

When the LCD (Liquid Crystal Display) low battery mark illuminates, the battery has become too weak to power the unit, and the existing battery should be replaced with a new battery. At this point, the unit will turn OFF until a new battery is inserted.

▲ NOTICE: GF Health Products, Inc. recommends the use of only a 9V alkaline or nickel-cadmium rechargeable battery with this device.

Replacing the Battery

When the LCD low battery mark illuminates, and the unit does not remain illuminated once turned on, the battery should be replaced.

1. Turn unit OFF.
2. Remove the front panel cover by pressing on the top of the panel and pressing down in order to slide the panel down. Continue sliding the panel downwards until the panel is completely removed from the unit. This will reveal the battery compartment.
3. Remove the discharged battery from the device.
4. Place new battery in the compartment. Note: Be sure the proper polarity (+ and -) markings match the markings in the device.
5. Replace the front panel cover.
6. Dispose of the old battery according to local guidelines and regulations.

CARING FOR YOUR DEVICE

Your device may be cleaned by wiping gently with a damp cloth moistened with mild soap and water. Do not immerse the device in water or other liquids.

Wipe lead wires with damp cloth moistened with soap and water. Do not immerse the lead wires.

To properly store the device for an extended period of time, remove the battery from the unit. Place the unit and accessories in the carrying case provided and store in a cool, dry location.

TROUBLESHOOTING

If the device does not function properly:

1. Ensure the battery is properly installed or replace the battery. Be sure to observe proper polarity markings when replacing the battery. If the LCD low battery mark illuminates when the unit is turned on, replace the battery and check again.
2. If the intensity has been adjusted and no stimulation is felt, check to ensure the lead wires are properly connected and the electrodes are properly applied to the skin. If the unit appears to be functioning and no stimulation is felt, the lead wires or electrodes may need to be replaced.
3. If the battery appears to be charged and the unit is not functioning, turn both Intensity Control Knobs to the OFF position (counterclockwise). Then gradually turn the Intensity Control Knob (clockwise) until stimulation is felt. If device still is not working, turn the unit off and contact your authorized GF Health Products, Inc. distributor.


If there is any other problem, please contact an authorized GF Health Products, Inc. distributor. Do not try to repair a defective device.

TECHNICAL SPECIFICATIONS

Channel	Dual, isolated between channels	
Pulse Intensity	Adjustable 0 mA – 80 mA peak into 500 ohm load each channel, constant current	
Pulse Width	30uS-300uS (adjustable)	
Pulse Rate	1Hz-150Hz (adjustable)	
Timer	10, 20, 30, 40, 50, 60 minutes and continuous	
Patient Compliance Counter	To show the usage time by user, counted per each 30 minutes treatment time	
Patient Lock System	To prevent changes in any set parameter by physician / clinician. To unlock or lock the unit, press the UP and DOWN buttons simultaneously for 2-3 seconds	
Wave Form	Mono positive-phasic square pulse, negative-phasic square pulse, asymmetrical bi-phasic square pulse	
LCD	Displays Timer, Function Mode, Alternate Mode, Pulse Width, Pulse Rate, Low Battery, Patient Compliance Counter, and Patient Lock System	
Function Mode	M Mode	Modulated width: Pulse width is automatically varied in an interval of 6 seconds. The modulation range of pulse width is from setting value to 60% less than the control setting value, then returns to the setting value. Pulse rate and pulse width are adjustable.
	MM Mode	Modulated width and rate: This mode consists of alternating modulated width and modulated rate. One parameter is always decreasing while the other is increasing and vice-versa. The modulation range is 50% decreasing.

Function Mode, Continued	MR Mode	Modulated Rate: The pulse rate is automatically varied in an interval of 0.5 seconds. The modulation range of pulse rate is from setting value to 50% less than the control setting value, then returns to the setting value. Pulse rate is adjustable but pulse width is fixed.
	B (Burst Mode)	Bursts consist of pulse width (adjustable), frequency=100Hz. Bursts occur twice every second.
	© (Constant Mode)	Pulse width and pulse rate are adjustable.
	SD Mode	Strength duration and width / rate modulation are adjustable. SD Mode consists of alternating modulated pulse width and pulse rate. One parameter is increasing while the other parameter is decreasing and vice-versa. The variety value is 50%; full cycle time is 22 seconds.
	Alternate Mode 	Stimulation output will begin from channel 2, last for 12 seconds, then switch to channel 1 for 12 seconds. During this mode, the stimulation will continue to alternate between the two channels.
Voltage	0 – 100 volt (open circuit)	
Power Source	9-Volt battery	
Battery Life	Approximately 70 hours at nominal settings	

LCD (Liquid Crystal Display)

The LCD displays timer (10, 20, 30, 40, 50, 60 minutes and © continuous), function modes (B, C, M, MM, MR, SD), alternate mode, pulse width, pulse rate, low battery, and patient compliance counter.

1. To check the LCD function, turn device ON and all the parameters will be displayed for 1-2 seconds.
2. After 1-2 seconds, the LCD will go to previous operation parameters.

Changing the LCD Parameters

- To set the treatment parameters on the TENS device, the unit must be in the UNLOCK position. If the device is in the LOCK position it, can be unlocked by pressing the UP and DOWN buttons simultaneously for 2-3 seconds.
- To switch between the different treatment parameters, press the MODE button.
- Once in the desired mode, press the UP, DOWN, or SET button until the desired treatment parameter is obtained.
- Once all of the desired treatment settings are displayed, the unit may be locked by pressing the UP and DOWN buttons simultaneously for 2-3 seconds while in the timer function.

LIMITED WARRANTY

SCOPE OF WARRANTY

GF Health Products, Inc. ("GF") warrants to the original purchaser only that it will replace or repair components, at GF's sole discretion, that are defective in material or workmanship under normal use and service. All warranties are conditioned upon the proper use of the products strictly in accordance with good commercial practice and applicable GF instructions and manuals, including proper use and maintenance. To the extent that a third party warrants a component, GF conveys all of its rights under that warranty to the original purchaser, to the extent permitted.

This limited warranty shall only apply to defects that are reported in accordance with the provisions set forth in this warranty document, within the applicable warranty period and which, upon examination by GF or its authorized representative, prove to be a warranty item. (See Obtaining Warranty Service below) This limited warranty is not transferable.

The warranted components and time periods are set forth below:

TENS Unit: 1 (one) year

* Labor is not included in the warranty.

† The warranty period is as designated above. If a part is replaced under warranty, the original warranty period will not be affected. All other replacement parts will be subject to the warranty period specified.

The applicable warranty period shall commence from date of shipment to the original customer, unless there is an expiration date on the component in which case the warranty shall expire on the earlier of warranty period or the expiration date.

OBTAINING WARRANTY SERVICE

Customers located in the United States who wish to report a warranty issue, must contact GF directly by calling 1.770.368.4700 or by e-mailing a request to cs@grahamfield.com. Customers located outside the United States must contact the Distributor from whom they purchased the products. In both cases, further directions will be provided once the initial contact is made. This limited warranty shall only apply to defects that are reported within the applicable warranty period. Failure to abide by the specific directions will result in denial of the warranty claim.

The warranty does not cover and GF shall not be liable for the following:

1. Defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering or failure to seek and obtain repair or replacement in a timely manner;
2. Products which are not installed, used, or properly cleaned and maintained as required in the official manual for the applicable product;
3. Products considered to be of a non-durable nature including, but not limited to: filters, fuses, gaskets, lubricants, and charts;
4. Accessories or parts not provided by GF;
5. Matching of color, grain or texture except to commercially acceptable standards;
6. Changes in color caused by natural or artificial light;
7. Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing, in advance, by GF;
8. Any labor or shipping charges incurred in the replacement part installation or repair;
9. Costs and expenses of regular maintenance and cleaning; and
10. Representations and warranties made by any person or entity other than GF.

ENTIRE WARRANTY, EXCLUSIVE REMEDY AND CONSEQUENTIAL DAMAGES DISCLAIMER

THIS WARRANTY IS GF'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. GF MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF ANY MODEL OR SAMPLE WAS SHOWN TO THE CUSTOMER, SUCH MODEL OR SAMPLE WAS USED MERELY TO ILLUSTRATE THE GENERAL TYPE AND QUALITY OF THE PRODUCT AND NOT TO REPRESENT THAT THE PRODUCT WOULD NECESSARILY CONFORM TO THE MODEL OR SAMPLE IN ALL RESPECTS. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF THE DEFECTIVE PARTS. GF SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO: DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, DOWNTIME, COVER, OR EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. CERTAIN STATES MAY CONFER ADDITIONAL RIGHTS REGARDING WARRANTIES AND IN THOSE STATES GF'S LIABILITY AND THE LIABILITY OF GF'S SUPPLIERS, SHALL BE LIMITED TO THE FULLEST EXTENT PERMITTED BY LAW.

The warranties contained herein, together with GF's current Terms and Conditions, contain all the representations and warranties with respect to the subject matter of this document, and supersede all prior negotiations, agreements and understandings with respect thereto. The recipient of this document hereby acknowledges and represents that it has not relied on any representation, assertion, guarantee, warranty, collateral contract or other assurance, except those set out in this document. For additional information on this product or this warranty, please contact a GF Customer Service Representative.

NOTES:

1. Additional terms and conditions may apply. See GF's General Terms and Conditions on its website: www.grahamfield.com.
2. Freight claims must be noted on the appropriate shipping documents and must be made with immediacy. International, federal and state regulations govern specific requirements for freight claims. Failure to abide by those regulations may result in a denial of the freight claim. GF will assist you in filing the freight claim.
3. Claims for any short shipment must be made within three (3) days of the invoice date.

Manufactured For
GF Health Products, Inc.
One Graham-Field Way, Atlanta, GA 30340

Made in Taiwan
© 2009, GF Health Products, Inc. All Rights Reserved.
GF is a registered trademark of GF Health Products, Inc.

GF_TXS-INS-LAB-Rev04

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website:
www.grahamfield.com

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