



### **Mounting the Bolsters:**

First remove tray. Next mount body bolster. Then mount head bolster.

Grip back bolster strap near attachment to cushion and stretch elastic strap. Slide into place over back cushion.

Position head bolster in the same way.

### **Removing the Bolsters:**

First remove tray. Next remove head bolster. Then remove body bolster.

Grip head bolster strap near attachment to cushion and stretch elastic strap. Slide up and off back cushion.

Remove back cushion in the same way.



### **Graham-Field Health Products**

2935 Northeast Parkway, Atlanta, Georgia 30360

telephone: 800-347-5678, 770-447-1609

fax: 800-726-0601, 678-291-3232

<http://www.grahamfield.com>