

## **Bariatric: Imperial Offset Cane**

- Heavy duty steel construction with 600 lb maximum weight capacity.
- Offset handle places the user's weight directly over the base for maximum balance and control.
- Attractive black finish resists marring and scratching.
- Easy push-button height adjustment with locking collet for added security.
- Height adjusts from 29" to 38" in 1" increments.
- Attractive black finish.
- Designed to meet the requirements of HCPCS code: E0100.
- Limited lifetime warranty.

6334A            4/cs



**6334A  
Black Finish**



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**[www.grahamfield.com](http://www.grahamfield.com)**