

Everyday Divided Leg Bathing Sling

- For use with four-point spreader option.
- Non-padded.
- Quick draining/drying polyester mesh make this sling ideal for bathing.
- Designed to meet the requirements of HCPCS Code: E0621.
- 6 month limited warranty.


400 lb. Weight Capacity

3303100	Small	1 ea
3303200	Medium	1 ea
3303300	Large	1 ea

250 lb. Weight Capacity

3301100	Small	1 ea
3301200	Medium	1 ea
3301300	Large	1 ea

Item	Sling Size	Body Width	Body Length	Total Width	Total Length
3303100, 3301100	Small	27 1/4"	14 1/4"	39"	50"
3303200, 3301200	Medium	37"	21 1/4"	57"	61"
3303300, 3301300	Large	40 1/4"	23"	60"	65"



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com