

Standard Ophthalmoscope Set

- Ophthalmoscope head with dioptic range of +20 to -20.
- Pinhole opening has light in diopetre scale.
- Chrome-plated screw lock handle with rheostat houses two C batteries (not included).
- Comes with replacement bulb and case.

1227-2 1/st



Graham-Field Health Products
2935 Northeast Parkway
Atlanta, Georgia 30360
2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com