

NEW AND IMPROVED!

Lumex Deluxe Three Position Recliner (Model 574RG)

Lumex continues to build on its heritage of innovation and product enhancement. Newly designed and updated, the Lumex Deluxe Three Position Recliners features all of the features of the standard recliners plus-

- Upgraded Tente® Casters. Two five inch rear swivel casters with wheel locks and two five inch front fixed casters for enhanced maneuverability in confined spaces.
- Easily removed “pop-out” side panels enable easy access to recliner for cleaning and maintenance.
- Welded steel frame is constructed of round tubing resulting in improved structural integrity of the frame.
- Improved wider and softer armrests for resident comfort.
- Weight Bearing Capacity: 290 Pounds

- **Lumex Wrap-Around Warranty**

1. Welded Frame: Limited Lifetime
2. Upholstery: Limited One (1) Year
3. All Other Components: Limited Three (3) Year

- **Available Colors:**

1. Rosewood: 574RG863
2. Blue Ridge: 574RG427
3. Jade: 574RG857
4. Grey Taupe: 574RG858
5. Doe Skin: 574RG851
6. Pecan: 574RG406
7. Terra Cotta: 574RG407
8. Royal Blue: 574RG454
9. Dove Grey: 574RG413

- **Fire Rated Colors**

**(Meets California Technical Bulletin 133
Flammability Standards)**

1. Rosewood: FR574RG863
2. Blue Ridge: FR574RG427
3. Jade: FR574RG857



Standard Colors



ROSEWOOD (863)



JADE (857)



BLUE RIDGE (427)



DOE SKIN (851)



GRAY TAUPE (858)



PECAN (406)



TERRA COTTA (407)



ROYAL BLUE (454)



DOVE GRAY (413)

- Upholstery color may vary slightly from brochure.
- Lumex Ortho-Biotic Seating with an FR prefix designates that the seating meets CA133 Fire Code specifications.
- Please contact Graham-Field Customer Service at 800-347-5678 for details on specialty colors and COM orders.



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com