

BP CUFF SIZE CHART

<u>BP CUFF</u>	<u>LIMB RANGE (inches)</u>
NEWBORN	3 <sup>3</sup> / <sub>4</sub> - 5 <sup>1</sup> / <sub>4</sub>
INFANT	5 - 7 <sup>1</sup> / <sub>2</sub>
CHILD	7 - 10 <sup>1</sup> / <sub>2</sub>
ADULT	9 <sup>1</sup> / <sub>2</sub> - 15
LG ADULT	13 - 19
THIGH	17 - 26