

- Virtually silent operation
- Mobile - with its small size, quiet operation and optional battery pack, the unit can be used almost anywhere
- Lightweight – compressor unit weighs only 0.55 lbs
- Flexible - can be powered using the AC power cord in a wall outlet, or the optional AC car lighter adapter or the optional battery pack
- 3 year limited warranty
- Designed to meet the requirements of HCPCS Code: E0574

JB0112-066	Nano-Sonic Nebulizer	1 ea
JB0112-066B	Nano-Sonic Nebulizer with Battery Pack, Auto Adaptor	1 ea
JB0112-065BATTERY	Optional Battery	1 ea
ADAP3011	Auto Adapter Only	1 ea

The John Bunn Nano-Sonic™ comes complete with:

- 1 pediatric mask
- 1 adult mask
- 1 AC adapter
- 1 mouthpiece
- 1 nasal piece
- 2 spare filters
- Optional portable, rechargeable battery pack
- Optional AC car lighter adapter



Technical Specifications

Method of nebulization	Ultrasonic vibration
Operating voltage	110V~ 120V~ 220V~ 230V~ 240V~ 50Hz 60Hz
Class	Class II
Power	18 W
Oscillation frequency	2.5 Mhz
Nebulization rate.....	0.5 ml./ min.
Particle dimension	1.5 - 5.7 (MMAD 3,8 micrones)
Medication cup maximum capacity.....	8 ml.

Weight	0.57 lbs / 260 g
Device dimensions	3.2" x 6.6" x 2.3"
Operating temperature and humidity	41°F to 104°F, 30% to 95% RH 5° to 40° C, 30% to 95% RH
Storage temperature and humidity	33.8°F to 104°F, 10% to 95% RH 1° to 40° C, 10% to 95% RH
Atmospheric pressure range for transport and storage	500-1060 mb



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com