



Feather® Sterile MiniScalpel

Sterile, disposable scalpel features heat-treated stainless steel blade securely integrated with rigid, chemically resistant plastic handle. This slender plastic handle is suitable for surgery in limited areas; such as brain, nerves, ears, nose, mouth, plastic surgery, blood vessels, skin and infants. Individually wrapped. Blades are protected with snap-off cover. Scalpels available with choice of 3 blade sizes.

All sizes packed 20/bx.

ISO 14001 compliant, ISO 9001 certified and ISO 13485 certified.

| | |
|---------|-------|
| 2978#11 | 20/bx |
| 2978#14 | 20/bx |
| 2978#15 | 20/bx |



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700019RevA07