

Digital Pediatric Scale

High performance digital pediatric scale with motion-sensing technology to capture the weights of even the most active babies.

- Weighs in kilograms and pounds/ounces
- Waterproof On/Off, Zero/Clear (tare), Hold/Release and Kg/lb buttons
- Features include auto-off, auto-zero and low-battery indicator.
- Tare memory
- Measuring tape included (24" x 1/8" / 59cm x 1mm)
- Powered by 6 AA batteries (included) or AC adapter (not included)

HOM553 1 ea



Display Panel



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700039RevA07