

Portamate™ II Oxygen Conserver

- State-of-the-art, light-weight, oxygen conserver.
- Features every-other-breath setting for up to 6/1 savings.
- Offers a continuous-flow toggle switch, pulse indicator light, low battery indicator, and built-in brass regulator with cylinder contents gauge.
- Operates on one C-size battery for 45 days of use.
- Built-in 0 – 3000 psi gauge.
- Three year limited warranty.

JB0170-001 1 ea

Width	Depth	Height	Weight	Power Requirements
2.8"	4.5"	1.5"	14 oz	One C-size Battery



Rear View



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700090RevA07