

**Lumex® Mark I Lift Chair**

Lumex®, the innovator in specialty healthcare seating introduces a brand new healthcare seating product designed for both the homecare and long term care markets. The Lumex Mark I Lift Chair has been designed for use in a variety of environments, including nursing facilities, physician offices, post-operative recovery rooms and the home. The Lumex Mark I Lift Chair functions both as a lift chair and as an adjustable height hip chair in a space that is about a quarter of the space that a typical lift chair would require.

- As a lift chair, the Lumex Mark I Lift Chair mechanism lifts the user from a seated to a standing position in a more vertical posture than standard lift chairs. This vertical motion helps the user maintain a more balanced posture when sitting or standing up from the chair.
- The direct drive motor with worm gear actuator allows for a very smooth and gentle lifting motion that is easily controlled by the ergonomic switch conveniently located under the front of the arm.
- As a hip chair, the reduced seat depth and the adjustable seat height of the Mark I facilitates the proper positioning of the patients body reducing stress on the hip and knee joints, crucial for proper healing. Unlike traditional hip chairs with stationary seat heights, the Lumex Mark I Lift Chair can be infinitely adjusted in seat height to properly fit the user.
- The reduced overall dimensions of the Lumex Mark I make it the perfect size for apartments, next to a bed or in a physician’s office. Two rear leg mounted wheels allow the Lumex Mark I to be easily moved within a room, between rooms or within a common area.

**GF4403-112** Rosewood Vinyl Upholstery 1 ea

**Options**

**Footstool:**

Convenient accessory to rest the users unaffected leg on when the Lumex Mark I is used as a hip chair. Can be stored under chair when not in use. Footstool is 6.25" in height.

**GF4403-FS12** Rosewood Vinyl Upholstery 1 ea

**Headrest:**

Comfortably supports users head. Can be easily repositioned with Velcro strip.

**GF4403-HR12** Rosewood Vinyl Upholstery 1 ea



**Specifications and Dimensions**

Minimum Seat Height:	19"
Maximum Seat Height:	27"
Seat Width:	18.75"
Seat Depth:	19.5"
Width Between Arms:	19.5"
Arm Height to Floor:	27"
Arm Height From Seat:	9"
Back Height From Seat:	28"
Overall Width:	24.5"
Overall Depth:	24"
Overall Height:	47"
Product Weight:	56 lb
Upholstery:	Complies CA117 Regulations
Power Unit:	IN:110 VAC OUT: 24VDC
Control Unit and Actuator:	UL/ CSA Approved
Weight bearing capacity:	264 pounds



**Graham-Field Health Products**

2935 Northeast Parkway  
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**