


Tissue Interface Pressure Map: Gold Care



mmHg	Distribution
50	5
47	0
44	0
41	1
38	0
35	1
32	1
29	0
26	1
23	3
20	6
17	29
14	34
11	37
8	44
5	78
2	2
0	2062

202 Cells = 5mmHg

mmHg	Distribution
50	0
47	0
44	0
41	9
38	29
35	36
32	46
29	45
26	62
23	65
20	97
17	95
14	112
11	98
8	87
5	111
2	1
0	1411

851 Cells = 5mmHg

mmHg	Distribution
50	3
47	1
44	2
41	1
38	3
35	8
32	10
29	40
26	78
23	81
20	93
17	89
14	70
11	97
8	106
5	138
2	4
0	1480

765 Cells = 5mmHg

General Comments
Gold Mattress June 10, 1999
Male 5' 6" 165 lbs