

Metro IC4

Metro IC4

The Metro IC4 is the all around lightweight wheelchair to fit both your product line and your budget. With over-link cross braces, flip back arms, height adjustable back, and anti-tippers, this high strength lightweight is the chair of choice. A lightweight wheelchair for a lightweight budget.



Metro IC4 (Industry Compatible)		Dimensions and Weights							
		Seat			Back	Overall			
			-lip	**			**		<u>⊥+</u> ∑
3D020120	Flip Back Desk Arm, Swingaway Footrest	18"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020130	Flip Back Desk Arm, Elevating Legrest	18"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020140	Flip Back Full Arm, Swingaway Footrest	18"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020150	Flip Back Full Arm, Elevating Legrest	18"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020220	Flip Back Desk Arm, Swingaway Footrest	16"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020230	Flip Back Desk Arm, Elevating Legrest	16"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020320	Flip Back Desk Arm, Swingaway Footrest	20"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020330	Flip Back Desk Arm, Elevating Legrest	20"	16"	19 ¾"	18 ½"	45"	38"	27"	11"
3D020340	Flip Back Full Arm ,Swingaway Footrest	20"	16"	19 ¾"	18 ½"	45"	38"	27"	11"

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". **With axles in lower position.



Graham-Field Health Products 2935 Northeast Parkway Atlanta,Georgia 30360 2007© Graham-Field Health Products

For Placing Orders: 1-800-347-5678 Fax: 1 www.grahamfield.com