

Walker Tray

- Allows you to carry personal items from room to room.
- Sturdy clip-on brackets attach to most standard walkers
- Easily folds down when not in use.
- Made of durable plastic making it easy to clean
- Raised edges prevents spills
- Features two cup holders.
- Easy to install. No tools required!
- Size: 16" (W) x 12" (D)

603900A 1 ea



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700175RevA07