

X-Ray Illuminators

Features:

- Designed for desk viewing or wall mounting.
- Each utilizes 15W fluorescent tubes for even diffusion of light.
- Housings are stainless steel with powder coated finish.
- Features Film Grip at top of unit with plexiglass viewing panel.
- Each bank designed with its own on/off switch for low power consumption.
- All units come with three-wire grounded cords.

3792GF	One Bank	17.3" x 22.8" x 4.7" (44cm x 58cm x 12cm)
3793GF	Two Bank	31.5" x 22.8" x 4.7" (80cm x 58cm x 12cm)
3795GF	Three Bank	45.7" x 22.8" x 4.7" (116cm x 58cm x 12cm)



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2007© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com

GF0700157RevA07