

CPR Pocket Size Resuscitator

- Ready, easy to use, compact system
- Ergonomic shape with soft air cushion that conforms easily to the face and provides patient comfort.
- Low resistance one-way valve minimizes risk of cross contamination
- Disposable filter and one way valve allows device to be reusable.
- Oxygen inlet allows product to be used with bag or any standard oxygen equipment
- For use on adults, and can be inverted for use on pediatric patients.
- Elastic strap keeps mask in place
- Folds into pocket storage case complete with CPR instructions

Latex Free

| | | |
|-----------------|---------------------|----------|
| GF91101 | Resuscitator Kit | 1 ea |
| GF91101F | Replacement Filters | 12 ea/bx |
| GF91101V | Replacement Valves | 12 ea/bx |



Specifications:

| | |
|-----------------------------------|----------------------------------|
| Dimensions: | 5.1" x 4.1" x 1.9" |
| Weight: | .2lbs (CPR mask with carry case) |
| Operation temperature: | 64°F – 104°F |
| Storage temperature: | 86°F – 122°F |
| Storage relative humidity: | <85% |

