

Traveler Recliner®

The same built-in quality as the standard Traveler wheelchair with the addition of a full reclining back. The back features unlimited adjustments using easy, finger-trigger hinge and cable mechanisms. The Traveler Recliner comes with a higher back and removable full support headrest. Rear axles are mounted further back than on standard Traveler for added stability. Anti-tippers standard. Meets Medicare K0001 / E1226 (Recliner Back) / E0966 (Head Rest Extension) specifications for wheelchairs.

- Maximum weight capacity 250 lb.
- 17" seat depth.
- Chair comes with Deep Blue upholstery.
- Burgundy and Black upholstery available as options. Please see Price List/Order Form for all available options
- 2 year limited warranty on frames and crossbraces.



Detachable desk arms
with elevating legrests

Traveler Recliner®		Dimensions and Weights								
		Seat		Back	Overall				Wt.*	
51012130	Detachable Desk Arm, Elevating Legrest	18"	17"	19 3/4"	23 3/8"	36"	53 1/8"	26 1/4"	13"	55 lb
51012150	Detachable Full Arm, Elevating Legrest	18"	17"	19 3/4"	23 3/8"	36"	53 1/8"	26 1/4"	13"	55 lb
51012230	Detachable Desk Arm, Elevating Legrest	16"	17"	19 3/4"	23 3/8"	36"	53 1/8"	24 1/4"	11"	55 lb
51012250	Detachable Full Arm, Elevating Legrest	16"	17"	19 3/4"	23 3/8"	36"	53 1/8"	24 1/4"	11"	55 lb
51012330	Detachable Desk Arm, Elevating Legrest	20"	17"	19 3/4"	23 3/8"	36"	53 1/8"	28 1/4"	13"	55 lb
51012350	Detachable Full Arm, Elevating Legrest	20"	17"	19 3/4"	23 3/8"	36"	53 1/8"	28 1/4"	13"	55 lb

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4". *Without Front Rigging.



Graham-Field
2935 Northeast Parkway
Atlanta, Georgia 30360
© 2008 GF Health Products, Inc.

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com