

Lo-Max™ Bed

- Mattress deck dimensions: 35" width and 76" or 80" mattress deck lengths.
- Sleep surface adjustable from 8 1/2" to 16 1/4".
- 400 lb weight capacity.
- Full-Electric or fixed height.
- Durable powder-coated epoxy finish.



Full-Electric Lo-Max Bed Shown with Spring Mattress Deck and Stock Medium Oak Finish, Rectangular T-mold Head and Foot Boards.

- Available with spring or pan mattress deck.
- Head articulation of 75 degrees for assistance with respiratory ailments.
- Foot/knee articulation of 45 degrees.
- Accommodates any head or foot board options.
- UL 2601-1 and CAN/CSA C22.2 No. 601-1-M90 approved.
- 76" and 80" mattress deck lengths.

Options and Accessories

Head and Foot Boards*	
AB53300	Rectangular, Solar Oak Finish For Full-Electric Bed
AB48300	Rectangular, Medium Oak Finish For Full-Electric Bed
AB80300	Rectangular, Montana Walnut Fin. For Full-Electric Bed
AB53350	Rectangular, Solar Oak Finish, Embedded Staff Control** For Full-Electric Bed
AB48350	Rectangular, Medium Oak Finish, Embedded Staff Control** For Full-Electric Bed
AB80350	Rectangular, Montana Walnut Finish, Embedded Staff Control** For Full-Electric Bed
AB29300CCR	Center Crown, Windsor Mahogany Finish For Full-Electric Bed
AB29350CCR	Center Crown, Windsor Mahogany Finish, Embedded Staff Control** For Full-Electric Bed
AB02300CCR	Center Crown, Figured Mahogany Finish For Full-Electric Bed
AB02350CCR	Center Crown, Figured Mahogany Finish, Embedded Staff Control** For Full-Electric Bed

* All boards above have T-mold edges.
Other styles and finishes available via custom order.
**Order in conjunction with AL75000 Embedded Staff Control

76" LENGTH	80" LENGTH	Mattress Deck Surface	Bed Type	Units
AZ4804	AZ4904	Spring	Full-Electric	1 ea
AZ4814	AZ4914	Pan	Full-Electric	1 ea

* Bed, head and foot boards, and rails sold separately.

Accessories	
554-2001-920	Emergency Crank**
554-0639-000	Trapeze Adapter
556-0005-000	Traction Frame Adapter (2 Req)
554-0143-000	Trapeze Bar (Patient Helper Bar)
554-0190-000	Panel-Mounted I.V. Holder
554-0769-000	Corner Mattress Retainer Kit
554-0765-000	Rail Clip On Pendant Holder

** Applicable to electric beds only.

Factory-Installed Options	
999-0700-355	Pull Pin Lock-Out Pendant
554-0658-000	Edema Foot Ratchet
AL75000	Embedded Staff Control

Side Rails and Assist Bars	
ZA7020	Half Head Rails (Set)
AZ7350	Assist Bar (1 Each)

Wall Saver Bumpers	
AZ63800	Wall Saver Bumpers, 76" Bed
AZ63400	Wall Saver Bumpers, 80" Bed

