





Panascope™ Stethoscopes-Lightweight

- This lightweight aluminum stethoscope offers a dual head chestpiece with a chrome-plated rotating stem.
- The non-chill ring is color coordinated with 22" latex-free tubing.
- Adjustable binaurals are chrome-plated brass with an internal spring.
- Comes complete with additional diaphragm and eartips.
- Lifetime Warranty

500	w/adult chestpiece	Black	1 ea	
500BL	w/adult chestpiece	Blue	1 ea	
500GY	w/adult chestpiece	Grey	1 ea	
510	w/pediatric chestpiece	Black	1 ea	
510GY	w/pediatric chestpiece	Grey	1 ea	