



**Now available with more durable, high-impact plastic bed ends!  
Attractive new plastic bed-ends are better able to withstand the storage, handling, cleaning, and use associated with homecare beds.**

The revolutionary single-motor design homecare bed is now available in a full-electric version. The WhisperLite™ II Full-Electric Homecare Bed offers the same features as the WhisperLite II Semi-Electric Homecare Bed with the additional feature of electric height adjustment. This makes it easier for caregivers to raise and lower the bed for patient care and linen changes. The electric height adjustment also enables the patient to easily position the bed for maximum comfort. Designed to be one of the lightest beds on the market. Foot section weighs only 53 pounds!



The WhisperLite II Full-Electric Bed is easier to deliver than any other full electric bed on the market.

- Total assembled bed weighs only 161 pounds (bed and motor).
  - Powerful DC motor for effortless height adjustment.
  - Motor comes pre-installed on foot section for easy set-up and disassembly.
  - All functions fully-integrated into convenient handheld pendant.
  - Welded motor-guard protects motor and universal joint from damage during shipment and deliveries.
  - Available with Fiberboard Decorative Walnut Bed-ends or High-impact Plastic Bed-ends.
  - Designed to meet the requirements of applicable HCPCS codes.
  - Maximum patient weight 350 lbs.
  - 450 lb evenly distributed total weight capacity including patient, mattress, rails, and accessories.
  - Overall Bed Size: 87" L x 36" W.
- Height Adjustment: Minimum Height: 15", Maximum Height: 24" (at deck level).
  - Lift Range: Head - 60°, Knee - 30°.

**Warranty:**

- Two year warranty for defects in workmanship and materials of mechanical components, frame, and electronics.
- During the warranty period, defective items will be repaired or replaced at manufacturer's option at no charge.

**Shipped in three cartons:**

1. Foot section and both motors.
2. Head section and Hi-Lo Rod.
3. Bed End Set.

Fiberboard Decorative Walnut Bed-Ends	High Impact Plastic Bed-Ends	Description	HCPCS
WL0458	WL0458PL	Full-Electric Bed Only	E0297
WL0458-220V	WL0458PL-220V	Full-Electric Bed Only - 220 V - International Only	
WL0458-PKG	WL0458PL-PKG	Full-Electric Bed with 1633-Innerspring Mattress and Full Rails	E0265
WL0458-PKGHLF	WL0458PL-PKGHLF	Full-Electric Bed with 1633-Innerspring Mattress and Half Rails	E0265
WL0458-BFPKG	WL0458PL-BFPKG	Full-Electric Bed with 1633-Basic Care Foam Mattress and Full Rails	E0265
WL0458-BFPKGHLF	WL0458PL-BFPKGHLF	Full-Electric Bed with 1633-Basic Care Foam Mattress and Half Rails	E0265
WL0458-BFPKGNAP	WL0458PL-BFPKGNAP	Full-Electric Bed with 1633-Basic Care Foam Mattress and No-Gap Rail	E0265
WL0458-FFPKG	WL0458PL-FFPKG	Full-Electric Bed with 1633-Extra Firm Foam Mattress and Full Rails	E0265
WL0458-FFPKGHLF	WL0458PL-FFPKGHLF	Full-Electric Bed with 1633-Extra Firm Foam Mattress and Half Rails	E0265
WL0458-IFPKG	WL0458PL-IFPKG	Full-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and Full Rails	E0265
WL0458-IFPKGHLF	WL0458PL-IFPKGHLF	Full-Electric Bed with 1633-Innerspring Mattress with Foam Sleep Surface and Half Rails	E0265
WL0458-RFPKG	WL0458PL-RFPKG	Full-Electric Bed with 1633-Reversible Foam Mattress and Full Rails	E0265
WL0458-RFPKGHLF	WL0458PL-RFPKGHLF	Full-Electric Bed with 1633-Reversible Foam Mattress and Half Rails	E0265
WL0458-XFIPKGR	WL0458PL-XFIPKGR	Full-Electric Bed with 1633-Extra Firm Innerspring Mattress and Full Rails	E0265
WL0458-XFIPKGHLF	WL0458PL-XFIPKGHLF	Full-Electric Bed with 1633-Extra Firm Innerspring Mattress and Half Rails	E0265
WL0458-RPKG	WL0458PL-RPKG	Full-Electric Bed with No Mattress and Full Rails	E0266
WL0458-RPKGHLF	WL0458PL-RPKGHLF	Full-Electric Bed with No Mattress and Half Rails	E0266
WL0458-RPKG-220V	WL0458PL-RPKG-220V	Full-Electric Bed with No Mattress and Full Rails - 220 V - International Only	



**Graham-Field**  
2935 Northeast Parkway  
Atlanta, Georgia 30360  
© 2008 GF Health Products, Inc.

For Placing Orders:

**1-800-347-5678**

**Fax: 1-800-726-0601**

**www.grahamfield.com**

GF0700031RevE08