

Replacement Bulbs

900	Otoscope And Anoscope Bulb 2.5v, .35a 8-32 Unc 2a Thread 15 Hours Life Tl-1 1/2 Nickel Plated, Vacuum	6/bx
901	Ophthalmoscope And Illuminator 3.5v, .72a 15 Hours Life T-1 1/2 Lamp Nickel Plated, Xenon	6/bx
902	Otoscope And Illuminator 3.5v, .72a 48 Grooves Per Inch 20 Hours Life Tl-1 1/2 Nickel Plated, Xenon	6/bx
903	Otoscope And Illuminator 2.5v, .64a 20 Hours Life Tl-1 1/2 Lamp Nickel Plated, Xenon	6/bx
904	Examination Bulb 14.5v, 2.3a 100 Hours Life T-3 Lamp 10-32 Tpi Thread Nickel Plated, Halogen, Xenon	6/bx
905	Child Laryngoscope Bulb 2.5v, .28a 1/8-72 Ns3 Thread 20 Hours Life Tl-1 Lamp Nickel Plated, Vacuum	6/bx
906	Adult Laryngoscope Bulb 2.5v, .38a 8-32 Unc 2a Thread 20 Hours Life Tl-1 1/2 Lamp Nickel Plated, Vacuum	6/bx
907	Ophthalmoscope Bulb 3.5v, .72a 20 Hours Life T-1 1.2 Lamp Nickel Plated, Xenon	6/bx
908	Fiberoptic Laryngoscope 2.5v, .66a 10-32 Unf 2a Thread 20 Hours Life Tl-1 1/2 Lamp Nickel Plated, Xenon	6/bx
909	Examination Bulb 15.1v, 2.26a 10-32 Tpi Thread 100 Hours Life T-3 Lamp Nickel Plated, Halogen, Xenon	6/bx
910	Vaginal Specula Illuminator 6.0v, .71a 10-32 Unf 2a Thread 75 Hours Life Tl-1 1/2 Lamp Halogen, Xenon	6/bx
911	Vaginal Specula Illuminator 4.6v, .81a 12-28 Nf Thread 30 Hours Life Tl-1 1/2 Lamp Halogen, Xenon	6/bx



900



901



902



903



904



905




906



907



908



909



910



911

