


Sure-Lift™ Padded Easy Lifting Sling with Headrest

- Provides excellent patient comfort and security
- Headrest provides additional support
- Constructed of high-performance Cordura that is more than twice as resistant to abrasions, tears, scuffs, and punctures than standard nylon, polyester, and canvas slings
- Reinforced back for support
- Padding on legs and back
- Height-adjustable straps
- Easy to use
- Machine washable and autoclavable
- For use with 4-point spreader bar
- One year limited warranty.
- Designed to meet the requirements of HCPCS code: E0621

Washing Instructions: Made of Cordura® with polyester tapes, padded with polyester foam. Can be washed at 80 degrees C and autoclaved at 85 degrees C for 30 minutes.

- | | | |
|----------------|--------|------|
| 3126100 | Small | 1 ea |
| 3126200 | Medium | 1 ea |
| 3126300 | Large | 1 ea |



Item	Sling Size	Tag Color	A	B	C	Weight Capacity
3126100	Small	Red	19"	29"	23"	350 lbs
3126200	Medium	Yellow	21"	33"	24"	660 lbs
3126300	Large	Green	25"	39"	25"	660 lbs

