

## Easy Raiser Belt


- For use with LF2011 and LF2025 Quick Raiser
- 660 lb maximum weight capacity
- 1 year limited warranty
- Designed to meet the requirements of HCPCS code: E0621

**Washing Instructions:** Made of Cordura® with polyester tapes, padded with polyester foam. Can be washed at 80 degrees C and autoclaved at 85 degrees C for 30 minutes.

<b>3120001</b>	Medium	1 ea
<b>3120002</b>	Large	1 ea
<b>31145109</b>	Optional Knee Strap	1 ea



*Optional Knee Strap*



Item	Sling Size	Tag Color	A	B	Weight Capacity
<b>3120001</b>	Medium	Yellow	13 ½"	34"	660 lbs
<b>3120002</b>	Large	Green	13 ½"	41 ½"	660 lbs

