

Non-Rebreather Mask

Reservoir bag with check valve, one safety vent open,
7' tubing, elastic strap and adjustable nose clip.

With Safety Vent

BF64060 Adult 50/cs

BF64090 Pediatric 50/cs

Without Safety Vent

BF64059 Adult 50/cs



Graham-Field Health Products

2935 Northeast Parkway
Atlanta, Georgia 30360

2005© Graham-Field Health Products

For Placing Orders:

1-800-347-5678

Fax: 1-800-726-0601

www.grahamfield.com