

Metro Power III

The Metro Power III is a transportable power wheelchair that provides an affordable alternative to scooters and high end chairs for qualified end users. It is a more traditional chair with effortless maneuverability both indoors and out.



- Rear wheels and casters with airless inserts
- Easy- folding frame for convenient transportability
- Highly reliable and industry-proven electronics by PG Drives Technology
- Ultra-durable frame that meets all ANSI-RESNA structural stress testing
- Available in black frame
- Value-added tension adjustable back upholstery for improved comfort and increased seating tolerance
- Convenient back pouch and front pouch for personal storage
- 300 lb weight capacity
- Warranty: Limited Lifetime on side frames and cross braces
- Medicare K0822 - Group II
- Batteries sold separately
- Additional models with batteries available (add "B" to the item number)

Metro Power™ III		Dimensions and Weights							
		Seat		Back	Overall			Wt.*	
2F100120	Desk Arms, Swingaway Footrest, Black Frame	18"	16"	19½"	17½"	41"	36½"	25"	84 lb
2F100130	Desk Arms, Elevating Legrest, Black Frame	18"	16"	19½"	17½"	41"	36½"	25"	84 lb
2F100220	Desk Arms, Swingaway Footrest, Black Frame	16"	16"	19½"	17½"	41"	36½"	24"	83 lb
2F100230	Desk Arms, Elevating Legrest, Black Frame	16"	16"	19½"	17½"	41"	36½"	24"	83 lb
2F100320	Desk Arms, Swingaway Footrest, Black Frame	20"	16"	19½"	17½"	41"	36½"	27"	85 lb
2F100330	Desk Arms, Elevating Legrest, Black Frame	20"	16"	19½"	17½"	41"	36½"	27"	85 lb

NOTE: Everest & Jennings maintains the policy of continual product improvement, therefore, we reserve the right to make changes without notice. All dimensions are ±1/4".

*Weight Capacity 300 lbs

*Without Front Rigging.

