


## Mesh Full-Body Commode Sling

- Quick draining/drying mesh fabric offers full head and neck support
- For use with 4-point spreader option
- Designed to meet the requirements of HCPCS code: E0621
- 6-month limited warranty

**Latex Free**

<b>DSLR114</b>	Medium	1 ea
<b>DSLR115</b>	Large	1 ea
<b>DSLR116</b>	X-Large	1 ea
<b>DSLR141</b>	Heavy Duty, XX-Large	1 ea



Item	Sling Size	Body Width	Body Length	Commode Opening to Top of Back	Total Length
DSLR114	Medium	42 3/4"	50 3/4"	34 3/4"	450 lb
DSLR115	Large	44 3/4"	54 1/4"	38"	450 lb
DSLR116	X Large	44 3/4"	60 1/4"	43"	450 lb
DSLR141	XX Large	42"	51"	41 1/2"	600 lb

