

Solid Full-Body Sling

- Fully padded fabric that offers a soft yet supportive surface
- Offers full head and neck support
- For use with 4-point spreader option
- Designed to meet the requirements of HCPCS code: E0621
- 6-month limited warranty

Latex Free

DSL112 Medium 1 ea
DSL113 Large 1 ea
DSL117 X-Large 1 ea



Item	Sling Size	Body Width	Body Length	Weight Capacity
DSL112	Medium	42 3/4"	50 3/4"	450 lb
DSL113	Large	44 3/4"	54 1/4"	450 lb
DSL117	X Large	44 3/4"	60 1/4"	450 lb

