



Stainless Steel Stethoscope

- Stainless steel dual head chestpiece has a non-chill ring for patient comfort
- The internal spring binaurals ensure acoustic sensitivity
- 22" latex-free tubing keeps acoustic sound intact, even in the noisiest settings
- Comes complete with two pair of contoured eartips
- Lifetime Warranty

LAB-7100	Black, Adult	1 ea	
LAB-7100GY	Grey, Adult	1 ea	
LAB-7200	Black, Pediatric	1 ea	
LAB-7300	Black, Infant	1 ea	
LAB-7300GY	Grey, Infant	1 ea	



LAB-7100