

# Alpha II™ & Növus™

The Therapeutic Bed




## One-Step Three-Point Safety Locking System

One-Step secures the bed in three locations – with just one step on the pedal – no need to walk around the entire bed to lock individual casters. One-step also easily unlocks the bed for cleaning or moving.


Saves time and adds security so that the bed won't move.



Shown unlocked and locked



Foot section unlocked and locked



Head section unlocked and locked



**OneStep AW7174030A**

80" Grid Deck, Full Electric with One-Step Three-Point Safety Locking System, with Staff Control, Backlit Handset, and Underbed Light

**OneStep NW7174**

80" Grid Deck, Full Electric with One-Step Three-Point Safety Locking System

**Basic American Alpha II™ Specifications**

Additional Features Include:

- Electrostatic powder coated antimicrobial bed surface provides longer lasting, durable finish
- Hospital Grade Grounded Plug is Standard

U.S. Patents

- US 7, 421, 748 1
- Multiple Patents Pending

Technical Specifications:

- Mattress Deck Length..... 80"
  - 4" extender available
- Mattress Deck Width..... 35" optional
  - 42" extender
- Overall Bed Length..... 85.5" with bumper
- Overall Bed Width..... 39.75" with rails/assist devices
- Maximum back angle ..... 70°
- Maximum knee angle ..... 30°
- Maximum multiple height ..... 30"
- Minimum height ..... 8.9"
- Horizontal translation ..... less than 1"
- Safe Workload (including resident, support surface, head/foot boards, and accessories) ..... 500 lb

**Basic American Növus™ Specifications**

Additional Features Include:

- High Quality UL 60601, IEC601-2-38, CAN/CSA C22.2 No. 601.1 compliance meets the highest global safety standards for electrical medical equipment
- Electrostatic powder-coated antimicrobial bed surface provides a longer-lasting, durable finish
- Hospital Grade Grounded Plug is Standard
- Comes standard with Corner Mattress Retainers

U.S. Patents

- US 7,421,748 1
- Multiple Patents Pending

Technical Specifications:

- Mattress Deck Length..... 80"
  - 4" extender available
- Mattress Deck Width ..... 35"/optional 42" extender
- Overall Bed Length ..... 85.5" with bumper
- Overall Bed Width ..... 39.9" with rails/assist devices
- Maximum back angle..... 70°
- Maximum knee angle..... 30°
- Maximum height ..... 30"
- Minimum height ..... 8.95"
- Horizontal translation ..... Less than 1"
- Safe Workload (including resident, support surface, head/foot boards and accessories) .....500 lb

For more information on the Alpha II or Növus beds, custom options and available accessories, please visit us online at [www.grahamfield.com](http://www.grahamfield.com) or call Basic American today.



**BASIC AMERICAN**  
MEDICAL PRODUCTS

770.678.4700

Specifications are subject to change without notice.  
Basic American Medical Products is a registered trademark, Alpha II and Novus are trademarks of GF Health Products, Inc. © 20011 GF Health Products, Inc. All rights reserved. GF1000122RevA11