

Tendercloud® Static Air Mattress Overlay

- Features over 240 flotation cells and chambers that self-adjust for pressure reduction
- Tri-layer system adjusts to redistribute air volume based on patient size and position
- Double-sealed heavy gauge 0.5mm vinyl to prevent air loss
- Emergency air release valve deflates overlay quickly and safely
- Four heavy-duty corner straps secure the overlay to the mattress
- Easily cleaned with mild detergent
- Instructions are on a welded flap to prevent loss
- Calendar included on flap to indicate date of beginning usage
- Compatible with most air filling devices
- 300 lb maximum weight capacity
- Designed to meet the requirements of HCPCS code: E0197
- 90 day limited warranty

Boston Air Valve

- Remove top black cap to fill the mattress (one way valve lets air in, but does not let any air out)
- Remove red and black cap to let air out
- Air valve is on extended flap to make it easy to add air anytime, even if someone is lying on the mattress

AQTC8035	72" x 35"	4/cs
AQTP-170	Manual Hand Inflation Pump	1 ea
AQRK-3	Tendercloud Repair Kit	1 ea

AQTC8035



AQTP-170