

## Cervical Philadelphia Collar

- Very high immobilization yet still provides patient comfort
- Molded out of high quality "Ethafoam" to ensure comfortable support
- Anatomical shape, preformed chin, and opening for a tracheotomy
- Sizes: Medium, Large, and X-Large

	<u>Size</u>	<u>Neck Height</u>	<u>Neck Width</u>
<b>8607M</b>	Medium	up to 70 mm (2 3/4")	16"-21"
<b>8607L</b>	Large	70mm-100 mm (2 3/4"- 3 15/16")	16"-21"
<b>8607XL</b>	X-Large	over 100 mm (3 15/16")	16"-21"

