

## Walkabout Lite Four-Wheel Rollator

The Lumex Walkabout Lite Four-Wheel rollator offers a comfortable, lightweight and affordable solution to mobility. Designed with an aluminum frame, this lightweight rollator weighs 14.5 lbs and supports up to a 300 lb weight capacity. The Walkabout folds quickly and easily into a compact unit for storage and transport. Other features include:

- Ergonomic hand grips
- Easy to operate locking and “quick stop” loop brakes (pull up to stop; push down to lock)
- Padded seat for added comfort
- Removable padded straight backbar with easy release pull knobs
- Adjustable handle height to accommodate broad range of user heights
- 6" wheel for indoor and outdoor use
- Basket included that can be placed underneath seat or in front for carrying of personal items
- Designed to meet the requirements of HCPCS codes: E0143 & E1056
- Includes easy to read, retail-friendly user manual
- Lifetime Limited Warranty

And now with added rear reflectors that increase visibility at night!

<b>RJ4300AQ</b>	Aqua	1 ea
<b>RJ4300B</b>	Royal Blue	1 ea
<b>RJ4300CH</b>	Champagne	1 ea
<b>RJ4300G</b>	Teal Green	1 ea
<b>RJ4300K</b>	Black	1 ea
<b>RJ4300L</b>	Lavender	1 ea
<b>RJ4300P</b>	Pink	1 ea
<b>RJ4300R</b>	Burgundy	1 ea

### Specifications

Width Inside Handgrips .....	16"
Overall Width .....	24.5"
Overall Depth .....	27.5"
Seat Width .....	14"
Seat Depth .....	13.5"
Seat Height .....	23"
Caster Diameter .....	6"
Minimum Handle Height .....	32"
Maximum Handle Height .....	37"
Weight Capacity .....	300 lb
Product Weight .....	14.5 lbs

*Specifications are subject to change without notice*



*Champagne model shown*



*Color options*

**Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.**



GRAHAM-FIELD

[www.grahamfield.com](http://www.grahamfield.com)