

Latex Free

Luminescent Sphygmomanometer w/ Gauge Guard

- Luminescent gauge
- Measure pressure from 0 mmHg to 300 mmHg
- Gauge guard to alleviate dropping impact on gauge
- Deluxe inflation bulb with filter screen
- Deluxe air release valve
- Nylon cuff, black
- Latex free
- 10 year calibration warranty
- 3 year warranty

203 Adult 1 ea
203X Large Adult 1 ea



<u>Cuff Class</u>	<u>Limb Range (inches)</u>	<u>Bladder (inches)</u>	<u>Cuff Dimensions (inches)*</u>
Newborn	3.75" - 5.25"	1" x 3.55"	7"W x 2"H
Infant	5" - 7.5"	2.35" x 4.72"	10.5"W x 3"H
Child	7" - 10.5"	3.55" x 5.90"	13.5"W x 4"H
Adult	9.5" - 15"	4.72" x 8.66"	20.5"W x 5.5"H
Large Adult	13" - 19"	5.71" x 12.41"	25"W x 6.5"H
Thigh	17" - 26"	7.09" x 14.57"	29.5"W x 8.5"H

* These are subjective measurements and can be vary ± 0.25"