

Folding Cane

- Cane folds to an 11 inch length for convenient transport and storage
- Derby-style handle is contoured for comfort
- Soft Grip models provide cushioning for added comfort to reduce hand stress and fatigue
- Easy push-button height adjustment from 32.5" – 36.5" with locking collet for added security
- 250 lb maximum weight capacity
- Attractive black frame finish resists marring and scratching
- Designed to meet the requirements of HCPCS code: E0100
- Limited lifetime warranty

5950A-4	Black	Standard Grip	4ea/cs
5950BL-4	Blue	Soft Grip	4ea/cs
5950G-4	Green	Soft Grip	4ea/cs
5950P-4	Pink	Soft Grip	4ea/cs
5950T-4	Tan	Soft Grip	4ea/cs

