

## Standard Sphygmomanometer

- Blue nylon cuff
- Latex two (2) tube bladder
- Standard inflation bulb and end valves
- Pressure measurement from 20 mmHg to 300 mmHg
- One (1) year calibration warranty
- Zippered carrying bag is included

**Caution:** This product contains natural rubber latex, which may cause allergic reactions.

**165NBL**    Adult    1 ea



<b>Cuff Class</b>	<b>Limb Range (inches)</b>	<b>Bladder (inches)</b>	<b>Cuff Dimensions (inches)*</b>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

\* These are subjective measurements and can be vary  $\pm 0.25"$