

Superior™ Sphygmomanometer

- Cotton cuff with gauge holder
- Color coordinated gauge and cuff – blue or grey
- Latex two (2) tube bladder
- Standard inflation bulb and end valve
- Pressure measurement from 20 mmHg to 300 mmHg
- Ten (10) year calibration warranty
- Zippered carry bag is included

Caution: This product contains natural rubber latex, which may cause allergic reactions.

175 Adult, Blue 1 ea
175GY Adult, Grey 1 ea



<u>Cuff Class</u>	<u>Limb Range (inches)</u>	<u>Bladder (inches)</u>	<u>Cuff Dimensions (inches)*</u>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

* These are subjective measurements and can vary \pm 0.25"

