

Labstar™ Deluxe Sphygmomanometer

- Accumax™ nylon cuff with index & range markings
- Latex two (2) tube bladder
- Large knurled air release valve for precision deflation
- Deluxe end valve reduces dust build-up
- Color coordinated gauge and cuff - blue, black, grey or orange
- Pressure measurement from 20 mmHg to 300 mmHg
- Lifetime calibration warranty
- Zippered carrying bag is included

Caution: This product contains natural rubber latex, which may cause allergic reactions.

202	Adult, Blue	1 ea
202C	Child, Blue	1 ea
202I	Infant, Blue	1 ea
202N	Newborn, Blue	1 ea
202T	Thigh, Blue	1 ea
202X	Large Adult, Blue	1 ea
202BK	Adult, Black	1 ea
202CBK	Child, Black	1 ea
202IBK	Infant, Black	1 ea
202TBK	Thigh, Black	1 ea
202XBK	Large Adult, Black	1 ea
202GY	Adult, Grey	1 ea
202XGY	Large Adult, Grey	1 ea
202COR	Child, Orange	1 ea
202IOR	Infant, Orange	1 ea
202OR	Adult, Orange	1 ea
202XOR	Large Adult, Orange	1 ea



202X



202GY



202BK



202



202T

<u>Cuff Class</u>	<u>Limb Range (inches)</u>	<u>Bladder (inches)</u>	<u>Cuff Dimensions (inches)*</u>
Newborn	3.75" - 5.25"	1" x 3.55"	7" W x 2" H
Infant	5" - 7.5"	2.35" x 4.72"	10.5" W x 3" H
Child	7" - 10.5"	3.55" x 5.90"	13.5" W x 4" H
Adult	9.5" - 15"	4.72" x 8.66"	20.5" W x 5.5" H
LG Adult	13" - 19"	5.71" x 12.41"	25" W x 6.5" H
Thigh	17" - 26"	7.09" x 14.57"	29.5" W x 8.5" H

* These are subjective measurements and can vary \pm 0.25"

